# add-i-gum

## Instructions for use / Please read carefully

### Composition:

Add-i-gum is a precision impression material based on addition-curing polyvinyl siloxanes, i.e. no volatile or gaseous by-products are created during setting.

**Putty:** ISO 4823:2015 type 0 – very high consistency – kneadable

**Light Body:** ISO 4823:2015 Type 3 – light bodied – low viscosity

### Shelf life:

Shelf life is 3 years, provided the material is stored at 23°C (74°F).  
Do not store over 25°C (78°F).

### Range of applications:

Add-i-gum can be used with both the two-step putty-wash and one-step putty-wash impression techniques. It is suitable for all areas of indication that require highly precise impressions.

### Additional information:

- Do not swap lids of the tubs or scoops of the putty as this could cause the material not to set.
- **Do not use latex gloves when handling the material!** The use of polyethylene or PVC gloves is recommended.
- Do not use with standard silicone impression materials (condensation-curing)!

### Pack contents:

#### Add-i-gum putty **535**

1 x 450 ml (680g) each of base (purple) and catalyst (white) plus 2 scoops

#### Add-i-gum Light N (Orange) **16306**

##### Normal setting

2 x 50 ml base and catalyst in safety cartridges,  
12 mixing tips and 12 intraoral tips

#### Add-i gum Light S (Green) **16339**

##### Fast setting

2 x 50 ml base and catalyst in safety cartridges,  
12 mixing tips and 12 intraoral tips

### Dispensing and mixing:

**Putty:** Remove equal quantities (according to volume or weight) of base (purple) and catalyst (white) paste from the tubs using the colour-coded scoops and knead until the material is homogeneous and has an even colour (45 seconds minimum).

Caution: close the tubs again immediately after use.

**Light Body:** Place the safety cartridge in the dispensing gun and attach a mixing tip. Syringe the preliminary impression with Light Body; for direct application to the tooth use intraoral tips. Take the wash impression without applying pressure.

### Technique:

Load the impression tray with the mixed material or apply the material directly to the tooth. The tray should be placed in the mouth no later than 1 1/2 min. (quick setting) or 2 1/4 min. (normal setting) after loading. Remove the tray after at least 2 min. (quick setting) or 3 1/2 min. (normal setting). These times apply for 23°C. Higher temperatures shorten the setting/working times while lower temperatures extend the setting/working times.

The material is thixotropic in the working stage, i.e. it only becomes flowable when pressure is applied during insertion of the tray or when it is extruded from the syringe. This allows the syringing of prepared teeth and enhances reproduction of detail in the subgingival region (to 5 µm).

### Disinfection

Impressions can be disinfected using standard disinfectant solutions.

### Pouring the impression

The impression should be stored at room temperature (max. 25°C).

Allow a minimum of 30 min. after removal from the mouth before pouring the impression. There are no further time restrictions.

### Electroforming

All impressions can be copper or silver electroformed.

### Technical data

#### Putty

#### ISO 4823:2015 type 0 – very high consistency – kneadable

Mixing time	45 seconds
Total working time (incl. mixing time)	2 minutes
Setting time (incl. mixing time)	3 1/2 minutes
Recommended intraoral setting time	3 minutes
Dimensional stability and final SHORE hardness	80 > 0.02%
Linear dimensional change after 24 h	< 0.1%
Elastic recovery	> 99.5%
Max. compression set	2.2%

#### Light Body

#### ISO 4823:2015 type 3 – light bodied – low viscosity

S=quick setting (green), N=normal setting (orange)

Total working time	1 1/2 min. (S), 2 1/4 min. (N)
Setting time	2 min. (S), 3 1/2 min. (N)
Recommended intraoral setting time	2 min. (S), 3 1/2 min. (N)
Dimensional stability and final SHORE hardness	50 > 0.05%
Linear dimensional change after 24 h	< 0.1%
Elastic recovery	> 99.7%
Max. compression set	6.5%

### Keep out of reach of children!

**Intended only for use by dentists.**

 **CB Healthcare Consulting GmbH**  
(manufacturer)  
**Am Neumarkt 34**  
**22041 Hamburg, Germany**

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