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CyberBond Adhesive LC

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CyberBond Adhesive LC

Single step bonding agent

CyberBond Adhesive is a dentin/enamel bonding agent which combines primer and adhesive into a single component. CyberBond Adhesive is compatible with all common dental etching materials.

INDICATIONS:

- Direct composite restorations.
- Adhesive cementation of crown & bridges (+ inlays + onlays).

ESSENTIAL INGREDIENTS:

- Multifunctional methacrylates
- Multifunctional monomers
- Urethan dimethacrylate
- Photoinitiators
- Acetone
- Ethanol
- Adhesive accelerator
- Colloidal Glass

INSTRUCTIONS FOR USE:

Preparation of the cavity:

Prepare the cavity of the tooth according to the clinical requirements.

Etching:

- 1) Etch enamel and/or dentin for 5-10 seconds. When bonding only to enamel is needed, etch enamel for 20-30 seconds. („CyberEtch“).
- 2) Rinse with a spray of oil-free water for 5-10 seconds.
- 3) Blot dry excess moisture from the surface. Verify that the surface has no visible moisture.

Application of CyberBond Adhesive:

- 1) Dispense directly onto a brush by holding the brush tip directly over the nozzle of the tilted bottle.
- 2) Apply ample amount of CyberBond Adhesive LC leave undisturbed for 15 sec. Remove solvents with a vigorous blast of air for at least 15 seconds or until no visible movement of the resin is noted. If the surface appears incompletely coated, apply additional material.
- 3) Cure CyberBond Adhesive LC for 20 seconds with a visible light-cure unit. (When the restoration retention area lies primarily within dentin, apply a second layer. Remove solvents immediately and light cure for 20 seconds).
- 4) Restore with visible light cured or self cure composite resin of your choice or cement restoration.

CAUTION!

CyberBond Adhesive LC contains polymerizable monomers which may cause allergic contact dermatitis in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.

STORAGE:

- Do not store above 24°C (75°F).
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- Tightly recap after use.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.

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