INSTRUCTIONS FOR USE

CyberCore • Radiopaque dual-cure composite for core build-ups

CyberCore is a dual-cure, fluoride releasing, resin composite material designed for the fabrication of core build-ups. CyberCore may be used with any Bis-GMA compatible bonding agent.

PROPERTIES:

- Easy to use
- Provides a consistent, homogeneous mix.
- Dual curing: material can be light cured, additional chemical cure assures that the material can be quickly bulk filled.
- Dense void-free core: cuts like dentin
- Fluoride release: strengthens natural tooth structure and helps prevent secondary caries.
- Radiopaque
- Superior compressive strength for long-lasting restorations.
- Enhanced mechanical properties of conventional Bis-GMA composite by utilization of intercalated and exfoliated nanocomposite technology

INSTRUCTIONS:

- 1) Prepare the tooth according to clinical requirements.
- 2) Isolate the tooth from oral fluids. A rubber dam is recommended.
- 3) Acid-etch the entire preparation with a 32-40 % phosphoric acid (like for instance CyberEtch) for 15 seconds.
- 4) Rinse off etchant with a water-spray for at least 15 seconds. Etched enamel should appear frosted.
- 5) Dry-etch dried surfaces.
- 6) Follow bonding agent instructions for use with dual-cure composites.
- 7) Assemble syringe mix tip. (Bleed syringe on initial use and bleed through assembled mixing tip and intraoral tip before each use).
- 8) Place intraoral tip of assembled syringe into base of prepared tooth and dispense an adequate amount of material to complete build-up. Let material rise around the tip while withdrawing to eliminate air bubbles.
- 9) Do not disturb during the gelation stage. Light cure for at least 40 seconds. Check for complete hardness. Chemical-cure intra-oral set-time is approx. 3 minutes from start to mix/application. Set-time takes longer outside of mouth.
- 10) Finish preparation with a bur.

TECHNICAL DATA:

Working Time:	1:30 minute	Solubility:	2 µg/mm3
Setting Time:	4:30 minute	Linear Shrinkage:	1.2 %
Compressive strength:	250 MPa	Radiopacity:	400 %Al
Diametral tensile strength:	40 MPa	Compatible with halogen light:	Yes
Flexural strength:	200 MPa	Compatible with plasma arc light:	Yes
Water sorption:	14 µg/mm3	Compatible with LED:	Yes

CAUTION!

CyberCore contains polymerizable monomers which may cause allergic contact dermatitis in susceptible patients. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.

PRECAUTIONS IN USE:

- Do not use in presence of eugenol based materials: complete setting of CyberCore may not take place.
- A calcium hydroxide or glass ionomer liner may be used to isolate CyberCore from eugenol sources.
- If refrigerated, allow to return to room temperature before use. Refrigeration will increase viscosity, making application more difficult, and will retard setting time.

ESSENTIAL INGREDIENTS:

- Triethylenglycoldimethacrylate • Multifunctional monomers • Polymerization activators • Photoinitiators • Glassfiller • Fumed Silica

STORAGE:

- Do not store above 24°C (75°F).
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.

According to the EU Medical Devices Regulation, users / patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

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CORE BUILD-UP AB
CB Healthcare Consulting GmbH