0473

## GB INSTRUCTIONS FOR USE

## DIRECTION &amp; SAFETY GB

## RECOMMENDATIONS FOR USE OF DEHP POLISHERS

All polishers are exactly tailored to their specific recommended usage. Improper use can lead to tissue damage, early wear or destruction of the polisher, but not mention the risks to the user, the patient or other persons.

## PROPER USE

These instruments must only be used with turbines, handpieces and contra-angles that are technically and hygienically perfect, meaning that they should be well maintained and correctly cleaned. The turbines and contra-angles used for these instruments must ensure precise and regular rotation. Insert the instrument shaft into the driving unit as far as possible. Before applying the instrument on the surface to be treated, bring it to the rotation speed needed. In order to avoid deformations, use circular movements with the instrument. Do not bend the instrument or use it as a lever because this would lead to increased risk of breaking the instrument. Immediately discard any deformed, non-rotating instruments. Unmounted polishers must be centered after shaftmounting in order to avoid vibrations during use and only high quality mandrels should be applied. Use mandrels of inferior or quality could break them and cause injuries. Protective glasses should be worn at all times. In case of improper use or material failure, the mandrel, shaft or workpiece can break and could become a dangerous flying object. Instead of wearing protective goggles, you can work behind a protective glass pane. Wearing a respiratory mask will prohibit inhaling dust from polishers or workpieces. We recommend a dust aspiration system.

IMPROPER USE OF THESE INSTRUMENTS LEADS TO BAD QUALITY RESULTS AND MORE RISKS. THIS IS WHY THEY MUST BE USED ONLY BY QUALIFIED PERSONS.

## INSTRUCTIONS REGARDING ROTATING SPEED

Never exceed the maximum rotation speed stated. The recommended and maximum rotation speeds allowed can vary between products. For that reason, please refer to the rotation speeds indicated on the packaging. In case of non-observance of the maximum rotation speed allowed, the polisher will tend to vibrate. Such vibrations can deform the polisher or its shaft and cause them to break. In the latter case, the user, the patient and other persons can be injured. Observing rotation speeds ensures optimum results.

## NON OBSERVANCE OF THE MAXIMUM ROTATION SPEED ALLOWS INCREASED SAFETY RISKS.

## PRESSURE ON THE POLISHER

Excessive pressure on the polisher can destroy it. Excessive pressure results in more heat produced. Excessive pressure can cause early wear of the polisher.

**ALWAYS AVOID EXCESSIVE PRESSURE ON THE INSTRUMENT BECAUSE THIS WILL CAUSE OVERHEATING, WHICH COULD DAMAGE THE PULP.** ADDITIONALLY, IN CASE OF VERY EXCESSIVE PRESSURE, THE INSTRUMENT CAN BREAK AND, AS A RESULT, PERSONS CAN BE INJURED.

## WATER COOLING

In order to avoid the overheating of a tooth, sufficient water cooling must be ensured (50 ml/min).

INSUFFICIENT WATER COOLING CAN LEAD TO IRREVERSIBLE DAMAGES TO THE TOOTH AND ITS SURROUNDING TISSUES.

## PREPARATION &amp; REPROCESSING GB

## WARNING NOTICES

Strong acids and strong bases may oxidize the stainless steel shaft. Rinse the polisher with distilled water after the treatment with de-cleaning and disinfectant solutions. Avoid temperatures > 130°C. Ultrasonic bath must not exceed temperatures of 42°C because of the possible coagulation of albumen.

## RESTRICTION OF REPROCESSING

Disposable products delivered unsterile, marked with the symbol may only go through the validated sterilization cycle ONCE before initial use.

## COMMENT

All instruments are delivered unsterile and must be run through the indicated cycle before and after each use. Be label on the cleaning and/or disinfecting solution must specifically say "suitable for rubber polishers or synthetics/silicones", because if the label only says "rotary instruments", the solution may not be suitable for polishers.

## INSTRUCTIONS

## CLEANING PREPARATION:

Pre-clean under running water with a brush (plastic) directly after use.

## CLEANING: MANUAL

Under running water with a brush (plastic bristles). Machine assisted: Ultrasonic assisted cleaning with a suitable cleaning agent and disinfectant. Rinse the instruments afterwards under running water.

## OR CLEANING: AUTOMATIC

Thermal disinfectant: Manufacturer's specification according to EN ISO 15883. Cleaning programme as indicated by the manufacturer in the operating instructions.

## DISINFECTION:

Use a solution classified as suitable for rubber and silicone polishers and synthetics by the disinfectant manufacturer. Exposure times and concentrations recommended by the manufacturer should be adhered to.

## RINSING:

Rinse the polishers with distilled water after treatment with cleaning and disinfectant solutions.

## DRYING:

Dry with fresh, clean, lint-free cellulose tissues.

DE Healthcare Ltd  
Unit 9, Kingsthorpe Business Centre,  
Studland Road, Kingsthorpe,  
Northampton NN2 6NE U.K.

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41 Weaver Road  
Denver, PA 17517 USA

## MAINTENANCE:

Visual check of all instruments with optical magnification (5-10 fold).

## INSPECTION:

No residues > continue to sterilisation.  
Visible residues > repeat cleaning.  
Reject and dispose of instruments in the event of discernible defects.

## STERILISATION:

For all instruments that need to be sterilized in accordance with EN ISO 17664 and all national valid legal requirements.

- Steam sterilization: Appliance according to EN 13060, validated procedure.
- C category: S or B-sterilizer.
- Holding time: Full cycle 5 min
- Sterilisation temperature: 134°C
- Drying time: 10 min
- Threshold values of contents for feed-water and steam condensates
- Loading of sterilizer according to manufacturer's instructions
- Follow manufacturer's operating instructions

## STORAGE:

Store instruments padded and protected from recontamination in proven suitable sterile packaging, cassettes or retainers.

## ADDITIONAL INFORMATION:

Repetitive reprocessing can change both the look and feel of the instruments slightly, but does not interfere with the instrument's function.

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