NANO-HYBRID COMPOSITE

• Instructions for use

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Contra-indications	<p>If the site cannot be isolated after enamel etching</p>
---------------------------	--

and during application and curing of DEHP Nano-Hybrid Composite.
If allergic to components of DEHP Nano-Hybrid Composite.
If oral hygiene is poor.

Side effects	<p>No systemic side effects are known. Contact allergies with products of similar composition have been reported in isolated cases.</p>
---------------------	---

Interactions with other agents	<p>Agents which contain eugenol and/or oil of cloves may affect the polymerization of DEHP Nano-Hybrid Composite. The use of zinc-oxide-eugenol cements in combination with DEHP Nano-Hybrid Composite should therefore be avoided. Discolorations may occur when using cationic mouth rinses, plaque indicators or chlorhexidine.</p>
---------------------------------------	--

Application	<p><i>Cleansing and isolation</i> Cleanse the tooth to be treated and the adjacent teeth. Select color and create a dry working area.</p>
--------------------	---

Preparing the cavity	<p>When preparing the cavity, every effort should be made to use a technique that is gentle to the tooth substance (principle of Adhesive Restoration Technique). Enamel edges in the anterior region are to be beveled. In the posterior region, only sharp enamel edges should be slightly beveled or rounded.</p>
-----------------------------	--

Pulp Protection	<p>A suitable pulp-protective liner should be applied in close proximity to the pulp. Procedure according to manufacturer's instructions for use.</p>
------------------------	---

Application of bonding agent	<p>Apply the bonding agent. Procedure according to manufacturer's instructions for use.</p>
-------------------------------------	---

Application of DEHP Nano-Hybrid Composite	<p>Dispense material out of the syringe or tip and apply into the cavity by using the increment technique.</p>
--	--

Shade	Thickness of the layer	≥800 mW/cm²	≥1600 mW/cm²	exposure time
DEHP Nano-Hybrid Composite A1/D2, A2/B2, A3/D3, A3.5/B3, A4/M5	2 mm	20 s	10 s	

Shade	Thickness of the layer	≥800 mW/cm²	≥1600 mW/cm²	exposure time
DEHP Nano-Hybrid Composite A1/D2, A2/B2, A3/D3, A3.5/B3, A4/M5	2 mm	20 s	10 s	

Maximum thickness of the layer and curing times:	<p>Maximum thickness of the layer and curing times:</p>
---	---

Shade	Thickness of the layer	≥800 mW/cm²	≥1600 mW/cm²	exposure time
DEHP Nano-Hybrid Composite A1/D2, A2/B2, A3/D3, A3.5/B3, A4/M5	2 mm	20 s	10 s	

Definition	<p>DEHP Nano-Hybrid Composite is a highly filled, radiopaque, universal nano-hybrid composite. DEHP Nano-Hybrid Composite covers a wide range of tooth shades with fewer shades.</p>
-------------------	--

Composition	<p>Methacrylates Barium glass, silanized Amorphous silica, hydrophobed</p>
--------------------	--

Technical Data	<p>Complies with ISO 4049 Range of dimensions of inorganic filler particles: 0,04–2,9 µm Inorganic filler content by volume: 58 % Inorganic filler content by weight: 77 %</p>
-----------------------	--

Indications	<p>DEHP Nano-Hybrid Composite is specifically indicated for:</p> <ul style="list-style-type: none">Direct filling of class I, II, III, IV and V cavities Luting and repairing of ceramic and composite restorations more translucent material in the layer technique
--------------------	---

Contra-indications	<p>If the site cannot be isolated after enamel etching</p>
---------------------------	--

Date of issue	<p>09-2017</p>
----------------------	----------------

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