# SPHERICAL AMALGAM

## • Non Gamma 2

Spherical amalgam

• INSTRUCTIONS FOR USE

Amalgama esférica

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Amalgame sphérique

• MODE D'EMPLOI

Kapselamalgam

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Amalgama sferica

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Amalgama esféricas

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**97-92560** 1 Spill yellow cap 400 mg alloy/capsule

**97-92561** 2 Spill orange cap 600 mg alloy/capsule

**97-92562** 3 Spill green cap 800 mg alloy/capsule

### FOR DENTAL USE ONLY

A usage dentaire uniquement

Solo para uso dental

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Apenas para uso profesional

**CAUTION: Product contains mercury.**  
**ATENCIÓN: Producto contiene mercurio.**  
**ATENCIÓN: La produit contient du mercure.**  
**ACHTUNG: Das Produkt enthält Quecksilber.**  
**ATTENZIONE: La prodotto contiene mercurio.**  
**ATENÇÃO: Produto mercúrio.**

**CB Healthcare Consulting GmbH**  
(manufacturer)  
 Am Neumarkt 34  
 22041 Hamburg, Germany  
 Tel: +49 40 656680  
63DE Rev.2.07.2019

### INDICATIONS FOR USE

Intended for the use as a dental restorative material in the treatment of dental caries.

#### Contra-indications for Use

Dental amalgams shall not be used:

- As retrograde filling material.
- In patients with a known history of allergy to mercury, identified by a patch test.
- In patients where the kidney has been weakened as a consequence of impaired kidney function.
- For dental treatment of deciduous teeth, of children under fifteen years and of pregnant or breastfeeding women, except when deemed strictly necessary by the dental practitioner based on the specific medical needs of the patient.

	1 Spill	2 Spill	3 Spill
Capsule Color*	Yellow Cap	Orange Cap	Green Cap
Identification of alloy content	400 mg	600 mg	800 mg
Mercury mass/capsule	290 mg	400 mg	530 mg

\*base color identifies the brand

Composition / Particle Shape	Ag	Sn	Cu	Zn	Hg	Spherical	Lathe Cut
Spherical Amalgam	56 %	28.8 %	15.2 %	-	40 %	100 %	-

- All percentages are by weight.
- The alloy to mercury ratio varies depending on the size and setting time.

#### Triturating Time:

Amalgamator	Speed	1 Spill	2 Spill	3 Spill
Silamat		4-5	5-6	7-8
Vari-Mix II	H	10-15	15-20	20-25
Duomat (4000rpm)	3200	5-6	5-6	7-8
ProMix	High (Rabbit)	5-7	7-9	9-12

- CyberFill amalgam gives satisfactory results with modern high speed vibrators.
- Not designed to be used with planetary mixers such as Copernicus and Capmix.
- Slow speed amalgamators should not be used

	Condense	Carving
Working time	2-3.5 min.	3.5-5 min.

- working time may vary, depending on triturating conditions and technique applied
- Polish 24 hours after placement.

#### Mechanical properties

Compressive Strength (24 hr): 375 - 475 Mpa

Dimensional Change: ≤ 0.15 % / Creep: ≤ 0.75 %

**NOTE:** The triturated amalgam should be bright and have a plastic consistency. If the amalgam is excessively plashy and wet looking and/or sets fast, decrease trituration time or amalgamator speed. A mix that appears excessively dry indicates that trituration time or amalgamator speed should be increased. Your clinical experience in obtaining a satisfactory working consistency will be the best guide to trituration times.

#### Self Activating capsules – Instructions for Use

- Place and secure the capsule in the mixing arms of the desired amalgamator.
- To triturate, set the desired speed and time of the amalgamator.
- After trituration, remove the capsule from the amalgamator. Before opening, capsule to release amalgamated mass from inner walls of capsule.
- Before opening, tap amalgam into the body of the capsule.
- With gloved hands, grasp the cap of the capsule with one forefinger and thumb, grasp the body of the capsule with the other forefinger and thumb, and pinch the cap while pulling apart. The cap should click apart.
- The amalgam in the body of the capsule should be tapped into an amalgam well.
- Inspect amalgam mass to ensure that the membrane from the mercury pillow pack is not mixed in with the amalgam. Remove the membrane from the amalgam if this is observed.
- Begin placement of amalgam.
- After placing the amalgam, snap the capsule body and cap back together and dispose of in accordance with all applicable regulations. Adherence to the ADA's current "Best Management Practices For Amalgam Waste (BMPs)" and "Dental Mercury Hygiene Recommendations" is strongly recommended.

#### Placement of amalgam

- Insertion of the amalgam should commence immediately after trituration.
- Vertical condensing should begin immediately after triturating using instrument with fiat smooth surface.
- Carving of the dental amalgam may begin immediately after completion of condensation.
- Relieve marginal ridge area with an explorer before removal of matrix band.
- A final polish will enhance the restoration's qualities. Conventional finishing and polishing techniques may be used.

#### Special Notes

- CyberFill alloys require less mercury than other alloys therefore, when triturating with dentomat adjust the regulating screw for a leaner mixture.
- Spillages: Mercury presents a health hazard if incorrectly handled. Spillages of mercury should be removed immediately according to local regulations.
- Moisture contamination: If moisture is introduced into the amalgam before it has set, properties such as strength and corrosion resistance may be affected adversely. If the alloy contains zinc, such contamination may result in an excessive expansion (delayed expansion). Whenever it is possible, use a dry field.

**Storage:** It is recommended that this product be stored in a well ventilated place.

**Expiration Date:** Do not use this product after its expiration date

#### WARNING - CONTAINS MERCURY

- The capsules and the dental amalgam product contain mercury.
- Mercury vapors may be harmful if inhaled.

#### Precautions

- The use of a rubber dam may decrease the amount of mercury absorbed by a patient during the removal or placement of an amalgam.
- Use this product with adequate ventilation (vapor exposure should be kept under 0.05 mg/m<sup>3</sup> or any lower limit set by any governmental or regulatory agency governing your use of the product).
- Do not place dental amalgams in direct contact with other types of metals.
- Dental amalgams are for single use only.
- Wear gloves when handling the self-activating capsules and dental amalgams.
- Use an enclosed amalgamator when mixing the self-activating capsules.
- Properly dispose of spent capsules and any excess unused amalgam.
- Use HGX, or similar type mercury-absorbing compounds in the event of spillage of the contents of the self-activating capsules.
- If a patient experiences a localized hypersensitivity reaction to dental amalgam, the amalgam should be removed.

### EN

These precautionary procedures should always be used in addition to procedures recommended by your local regulatory agency and dental association.

#### Adverse Reactions

Effects and Hazards of Eye Contact: Irritant. Acute Exposure: Contact may cause irritation. Mercury is corrosive and may cause corneal injury or burns. Chronic Exposure: Mercury may be deposited in the lens of the eye, causing visual disturbances. Effects and Hazards of Skin Contact: Irritant/Sensitizer/Neurotoxin/Nephrotoxin. Acute Exposure: May cause redness and irritation. Chronic Exposure: Possible sensitization, dermatitis, and swelling. Mercury may be absorbed through the skin causing urinary problems. Effects and Hazards of Inhalation: Irritant/Sensitizer/Neurotoxin. Acute Exposure: Inhalation of mercury vapor can cause cough, fever, nausea, and vomiting. Chronic Exposure: Inhalation of high concentrations of mercury vapor over a long period may cause mercurialism. Findings are extremely variable and include tremors, salivation, stomatitis, loosening of teeth, blue lines on gums, pain, and numbness in extremities. Effects and Hazards of Ingestion: Neurotoxic/Nephrotoxic. Acute Exposure: May cause nausea, vomiting, kidney damage and nerve effects. Chronic Exposure: Symptoms include central nervous system (CNS) disorders.

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