Bulk Fill Universal

Light-curing posterior filling material

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
DE Healthcare Products
DEDICATED TO EXCELLENCE

DIN EN ISO 4049

CB Healthcare Consulting GmbH
Am Neumarkt 34
22041 Hamburg, Germany
Rev 2020/08

Light-curing, posterior filling material In accordance with DIN EN ISO 4049

• INSTRUCTIONS FOR USE

Bulk Fill Universal is a light-curing, highly radiopaque hybrid composite designed specifically for time-saving costeffective use in the posterior region. The combination of a new multi-hybrid filler technology with an innovative initiator system formed the basis for a filler material exhibiting minimal polymerisation shrinkage and excellent depth of cure. Bulk Fill Universal can even be cured reliably in increments of 4 mm and with very short polymerisation times. **Bulk Fill Universal** undergoes reduced shrinkage stress to ensure restorations with integral margins. The physical strength and high abrasionresistance are high enough for occlusal load-bearing fillings. **Bulk Fill Universal** contains 86% by weight inorganic filler (= 70.1% by volume) in a methacrylate matrix (Bis-GMA, UDMA, TEGDMA). **Bulk Fill Universal** is cured with halogen or LED light with a blue spectrum and is available in practical syringes or caps for direct application

Shade:

Universal (U)

Indications:

- Class I and II posterior fillings
- Core build-up

Preparation:

Clean the tooth to be treated. Mark the occlusal contact areas - minimal separation simplifies creating proximal contact areas and placing the matrix.

Allow the material to reach room temperature before use.

Cavity preparation:

It is essential to prepare the tooth minimally-invasively (adhesive fillings), preserving healthy tooth structure. Clean and dry the cavity. Follow the special requirements for bonding agents (see the relevant instructions).

Creating a dry workfield:

Ensure that the work area is sufficiently dry. Use of a rubber dam is recommended.

Matrix placement:

Take appropriate measures for areas where the applied Bulk Fill Universal can flow away. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and placement of the matrix.

Pulp protection:

A suitable pulp protective liner that can be covered with a stable cement, if necessary, should be applied in close proximity to the pulp.

Bonding material:

Bulk Fill Universal is used in the adhesive technique with a dentine-enamel bond. Any light-curing bond material can be used. Follow the respective instructions for use with regard to preparation (etching technique) and application.

Application of Bulk Fill Universal:

Prior to use, allow the material to attain room temperature. Proximal cavities require a matrix. It is better to use translucent matrices (smooth surfaces). Apply the Bulk Fill Universal and adapt it with a suitable instrument.

Bulk Fill Universal Caps:

Insert Bulk Fill Universal Caps into the opening of the dispenser (follow instructions for use). Turn the cap into desired direction and remove the protective top. Apply Bulk Fill Universal directly into the cavity by slowly and evenly pressing the levers of the dispenser together.

Light-curing:

Bulk Fill Universal is specially formulated for applying in layers of up to 4 mm.

Higher fillings must be built up in layers and each layer cured separately with a polymerisation lamp:

10 s when using lamps rated at 800 mW/cm² or higher.

20 s when using lamps rated at 500 to 800 mW/cm².

The outlet of the light guide should be positioned as close as possible to the restoration. Ensure that every section of the filling is fully exposed to the light. An additional polymerisation of the accessible case from the oral or vestibular side is recommended. Remove metal matrices before hand. If the light outlet is more than 5 mm from the surface, it must be assumed that the depth of cure will be reduced and that only that area actually illuminated by the light cone will be cured. Inadequate curing may lead to discolouration and pulpal complications.

Finishing:

After removing the matrices, the filling can be trimmed and polished immediately while being cooled (e.g. fine or extra-fine grit diamonds, polisher). The tooth should be fluoridated as a final step.

Information/precaution:

No known side-effects. Sensitisation in hypersensitive patients, however, cannot be excluded.

- **Bulk Fill Universal** contains Bis-GMA, UDMA, BHT, TEGDMA). **Bulk Fill Universal** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

- Phenolic substances, especially preparations containing eugenol or thymol, interfere with the curing of filling composites.

Avoid the use of zinc oxide eugenol cements or other eugenol containing materials in combination with filling composites.

- Fillings exposed to occlusal-loading should be examined at least once per year to ensure that any changes can be identified in time.

- Allow the caps to attain room temperature before use and apply with even, but not excessive, pressure.

Storage:

Store at temperatures of 4° C to 23° C. May be refrigerated. Do not expose syringes to any light. Do not use after expiry date. Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/ or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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