# BA Powersurge

IMPLANTOLOGY CONTRA-ANGLE 20:1  
CONTRE-ANGLE D'IMPLANTOLOGIE 20:1  
CONTRAANGULO PARA IMPLANTOLOGIA 20:1  
CONTRANGOLO D'IMPLANTOLOGIA 20:1  
IMPLANTOLOGIE WINKELSTÜCK 20:1  
CONTRA-ÂNGULO PARA IMPLANTOLOGIA 20:1

**PS20X - PS20L - PS20LED**

INSTRUCTIONS FOR USE - NOTICE D'INSTRUCTION  
NOTA INFORMATIVA DE INSTRUCCIÓN - FOGLIO D'ISTRUZIONI  
GEBRAUCHSANWEISUNG - INSTRUÇÕES DE USO

**PS20X / PS20L / PS20LED**  
**BA660070 / BA660570 / BA660571**  
REF 9799548 / REF 9799549 / REF 9799550

**B.A. INTERNATIONAL** 

## TABLE OF CONTENTS

- I. MEANING OF SYMBOLS
- II. SCOPE OF USE
- III. GENERAL SAFETY INFORMATION
- IV. TECHNICAL CHARACTERISTICS
- V. INSTALLATION OF MD
- VI. HYGIENE AND MAINTENANCE
- VII. REPAIR
- VIII. GUARANTEES
- IX. ACCESSORIES
- X. DISPOSAL

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

Sterilisable up to 135°C

Visual inspection

Medical device  
(contra-angle, hand  
piece, etc.)

General information

Wear gloves

Fibre optic device for  
motor with a light source

Device reference

Device with light  
emitting diode

Device serial number

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Date of manufacture

Thermo-disinfectable

Name of manufacturer

Single use

## II - SCOPE OF USE

Surgical treatment for hard organic substances: drilling, tapping, screwing of dental implants, screwing of the prosthetic parts.

Inappropriate or indirect use could damage this device and be a risk to the user and third parties.

These **MDs** are for professional use only in the dental surgery field.

According to these instructions, the **MD** must only be used by a user with dental medicine experience, for the stated use, and in compliance with current guidelines regarding the prevention of accidents at work and work protection and the instructions in this leaflet. **MDs** must be prepared and maintained only by individuals who have been trained in preventing infections, auto-protection and patient protection.

According to these instructions, users must :

- Only use non-defective work instruments,
- Only use the MD in a unit that complies with the guidelines of the EN 60601 standard,
- Follow the correct use instructions,
- Protect oneself and patients and third parties against all dangers,
- Avoid all contamination by the product.

The following situations:

- Inappropriate use,
- Lack of maintenance,
- Use of removable accessories or parts not approved by B.A. International,
- Using accessories from other devices on this **MD**,
- Change or addition to a **MD** not validated by B.A. International.

Relieve B.A. International of all guarantee obligations or other claims.

These **MDs** comply with the Community Directive 93/42/CEE as amended by the 2007/47/CEE directive.

### III - GENERAL SAFETY INFORMATION

Prior to use, check that the device has not been damaged in any way and that no parts are missing.

Wear appropriate protection, especially gloves, a mask and glasses.

Do not use the **MD** and inform your distributor or B.A. International After-Sales (AS) in the following situations :

- visible failure or damage,
- non-existent or insufficient irrigation (at least 50 ml/min, correctly tempered and oriented in an optimum manner),
- scratch or damage to the emitting side of the LED.

Risk of damage to the device and injury, wait until the motor has completely stopped to:

- connect/disconnect the **MD** from the motor,
- turn the assembly/disassembly system on,
- turn the locking/unlocking system for the rotary instrument on,
- handle the rotary instrument.

Assess the risk of septic substances penetrating tissue if used with a central air cooling motor (refer to the motor manufacturer's leaflet).

Warning : this **MD** has a ball bearing system. The power is maintained :

- For motors having a torque calibration system (like the **Powersurge** unit): you must calibrate the contra-angle as directed by the manufacturer. The motor must have been calibrated for the used contra-angle.
- For any other motors: make sure that the displayed torque values are right.

EN

FR

ES

IT

DE

PO

A tapping or a screwing torque which is too high (above 35N.cm) may increase risks of necrosis or of bone fissuring. Furthermore, the implant and the instrument may be damaged: implant manufacturer's prescriptions have to be respected.

Risk of burning and injury :

- Apply pressure to the **MD** with a rotary instrument (refer to the manufacturer's instructions) without exceeding 2N (about 200g).

Risk of burning related to the LED overheating :

- Avoid all contact between the LED and soft tissue.

Risk of dazzle:

- Do not look at the light output directly (especially on the LED model).
- Do not point the beam directly at the eyes of the patient or third parties.

The lighting level can generally be adjusted in the dental unit.

Risk of electrical discharge due to **MD** being badly connected to a B.A. International non-compatible system:

In the event of the **MD** being assembled and used in other manufacturer's treatment devices and installations, refer to "Protection against the risk of electrocution", "Leaking current" and "Non-earthing of use part" according to the IEC 60601-1 standard.

Do not use in explosive atmospheres.

## IV.2 - Characteristics of MDs

References	PS20X	PS20L	PS20LED
<b>Ratio</b>		20 : 1	
<b>Colour code</b>		Green	
<b>Adjustable depth stop system</b>	No	No	No
<b>Light</b>	No	Yes (light guide)	Yes (LED in the contra- angle)
<b>Motor connection standard</b>	ISO 3964	ISO 3964 <sup>1</sup>	ISO 3964 <sup>2</sup>
<b>Maximal speed of motor (rpm)</b>		40 000	
<b>Maximum torque of motor (N.cm)</b>		6,5	
<b>Maximum torque of instrument (N.cm)</b>		80 <sup>3</sup>	
<b>Rotary instruments according to EN ISO 1797-1</b> -> Ø (mm) -> Lg max (mm) -> Ø max in the active part (mm)		Type 1 2,35 40 <sup>4</sup> 10 <sup>4</sup>	
<b>Spray water flow according to ISO 7785-2 (ml/min)</b>		> 50	

S3 10% intermittent mode (1 min on / 9 min off). The authorised relative load duration is 1 minute per 9 minutes of pause.

<sup>1</sup> INTRAmatic LUX® connection system - Kavo® registered trademark. • <sup>2</sup> Powersurge LED B.A. International connection system. E-PLUS or E-PLUS Intra specific connection. • <sup>3</sup> The torque can be reached providing the motor used is powerful enough: for example, the Powersurge motor by B.A. International. In general, it is possible to limit the motor torque on the operating unit. • <sup>4</sup> Indicative values. In the event of the use of longer or bigger rotary tools, the user is responsible for choosing the right operational conditions that avoids all risks to patients or third parties. Follow the instructions for use, in accordance with the tools manufacturer's instructions.

EN

FR

ES

IT

DE

PO

#### **IV.3 - Adaptor for LED MD**

Intended to be placed instead of the bulb in dental motors with a INTRAmatic LUX® 1 type connection, the adaptor for LED **MD** will be supplied with power from the treatment unit.

Electrical characteristics for power of adaptor and/or LED **MD**:

- Nominal supply voltage for the adaptor for LED **MD**: 3.4 V CC
- Nominal supply voltage range for the adaptor for LED **MD**: 3.0 – 3.6 V CC
- Maximum current: 400mA

The LED is a semi-conductive element which must only be powered with continuous voltage. In order to ensure adequate operation, it must be properly polarised: the positive terminal of the connection part for powering the LED at the back of the **MD** is closest to the centre (**Fig 3**).

### **V - INSTALLATION OF MD**

The **MD** is supplied non-sterile.

Before first use, the **MD** must be cleaned, decontaminated and sterilised (see § **VI**).

#### **V.1 - Installation of the adaptor for the LED MD**

Check that your treatment unit satisfies the electrical power for the supply voltage (see § **IV.3**). If the electrical power is not indicated, please contact the manufacturer of your treatment unit or your installer.

- Disconnect the electrical supply for the treatment unit when installing the adaptor for the LED **MD**. Risk of electric shock.
- Proceed as with replacing the bulb – refer to the motor manufacturer's instructions.
- Insert the adaptor up to the stop.

The polarity of the power for the LED is important: if this does not light up when in use, then use an adaptor for the LED **MD** with a polarity correction or contact your installer.

#### **V.2 - Connection to the motor**

Check that the **MD** is completely dry before connecting it to the motor. Install the **MD** on the motor connection until ratcheting. For this, keep the motor and the **MD** in the same axle.

or

Turn the **MD** until the retractable catch lodges in the E notch of the motor (**Fig 4**).

Apply light traction to the **MD** to check that it is properly attached to the motor connection.

Operational test by switching on the motor; start at low speed, then gradually increase. If you notice overheating, irregularities, vibrations, abnormal noises when using the **MD**, immediately contact your dealer or your B.A. repair center.

#### **V.3 - Connection of the irrigation (**Fig. 5**)**

Attach the spraying nozzle to the head of the **MD**.

Connect the extremity of the irrigation pipe to the nozzle entry.

#### **V.4 - Disconnection of motor**

Remove the **MD** by keeping it in the motor axle.

In the event of prolonged non use of the device, do not leave it connected to the motor. Risk of damaging the motor.

#### **V.5 - Connection/disconnection of rotary instrument**

For all handling of rotary instruments, wear protective gloves. Risk of cutting and contamination.

Only use rotary instruments in good condition. Risk of necrosis.

**EN**

**FR**

**ES**

**IT**

**DE**

**PO**

## VI - HYGIENE AND MAINTENANCE

Sterilisation of medical devices must be done by properly trained and protected staff, in compliance with current regulations. The sterilisation protocol must be appropriate to the infectious risk.

Wear appropriate protective clothing: risk of infection and injury.

Only use products for the maintenance of medico-surgical equipment compatible with stainless steel (no chlorinated content).

Forbid antiseptics that are intended for use on skin and mucous.

Forbid products containing aldehyde, alcohol or other products likely to bind proteins.

For each product used: refer to the manufacturer's instructions. Comply with the concentrations, exposure durations, and life span of products. Do not mix the products and follow the instructions for their disposal.

### VI.1 - Preparation of MD in the place of use

To be done immediately after surgery (maximum time of 2 hrs):

- Disconnect the **MD** from the motor,
- Remove the rotary tool (see § VI.6),
- Remove the largest organic contaminants with disinfectant wipes.

- Wipe a disinfectant wipe over each of the **MD**'s sub-units,
- Or Spray disinfectant on the **MD** and wipe with a clean cloth.

#### Or: by Ultrasound

- Use a low frequency ultrasound tank (25 to 50 kHz) and a detergent-disinfectant product compatible with this process and with the **MD**,
- Friction of parts between them or with the tank can cause appearance defects,
- Install light input/output protections on the body (**Fig 7**) to protect the fronts of the fibre optic or the emitting surface of the LED,
- Immerse the components.

#### Rinsing and drying

- Empty the spray pipe (**Fig 4**).
- Rinse well<sup>5</sup>, then dry each part.

#### VI.2.2 - Automatic preparation for sterilisation

- Only in a washing and disinfecting machine, appropriate for the treatment of this type of **MD**. Follow the manufacturer's instructions for the device.
- If the device is not equipped with pipes for contra angles / handpieces, then disassemble the **MD** (see § V.2) and immobilise each part. Otherwise, put the **MD** mounted on the pipe.
- The disinfection thermo cycle must be at least 10 minutes at 93°C (203°F).

- Check there is no residue and that all the parts of the **MD** are completely dry at the end of the cycle.

- If the washing and disinfecting machine is near the place of use, it is possible to also do step VI.1.

#### Re-assembly of MD:

Only assemble parts with the same serial number<sup>6</sup>.

- Attach the collar cartridge D to A head (only one direction possible) until it completely plates (turn the cog until connection if required),

- Insert the head – collar cartridge unit in the G sleeve (connect a clean tool in the head and turn it slightly so that the cog connects),
- Bring the F lever to the central position, until ratcheting. Make sure that the depth stop system is in the right position.

Systematically check the condition of the front part by applying axial traction.

EN

FR

ES

IT

DE

PO

13

<sup>5</sup> Temperature < 38°C, distilled water is recommended instead of running water if this has too much chlorine (cf. FD98-135 standard).

<sup>6</sup> Exception : If the After-Sales Service has to replace a part, the serial number may be different.

### VI.3 - Lubrication

Follow the spray lubricant manufacturer's instructions.

Wear appropriate protective clothing (splashes, etc.).  
Keep away from all sources of heat or ignition.  
In particular: do not smoke. Risk of inflammation.

After cleaning and before each sterilisation (once a day minimum):

- Remove the rotary instrument,
- Insert the tip to the back of the **MD** (**Fig 8**),
- Cover the head of the **MD** with a soft cloth or paper or a wipe.
- Point the head towards the bottom,
- Spray several times by firmly holding the instrument,
- Wipe excess oil with a cloth or a wipe.

### VI.4 - Operational test

- Connect the **MD** to a micro-motor, point the head towards the bottom.
- Turn the micro-motor on at low speed for 30s. Gradually bring up to maximum speed.
- Wipe excess oil with a cloth or a wipe

If you notice overheating, irregularities, vibrations, abnormal noises when using the **MD**, immediately contact your after-sales service.

### VI.5 - Sterilisation

- The instruments are to be sterilised before first use and after each use.
- Only sterilise cleaned, lubricated and tested instruments.
- We recommend sterilisation that complies with the EN 13060 standard, class B. All other methods are forbidden.
- Independently put each **MD** in a sterilisation sachet that complies with current standards and the size of the **MD**.
- Comply with the space between the sachets and do not overload the autoclave.

- These **MDs** and their accessories (apart from accessories with the logo which are single use) must be sterilised at 135°C for 18 minutes minimum (sterilisation time).

- Check that the device is completely dry at the end of the cycle.

### VI.6 - Storage

Keep the **MD** in a sterilisation sachet in a clean and dry place.

## VII - REPAIR

In the event of breakdown, please contact your approved distributor or your local B.A. International repair centre.

Repairs must only be carried out by an approved repairer or by B.A. International After Sales Department, only with B.A. International original replacement parts.

For all revisions or repairs, the **MD** must be returned complete and sterile with proof of sterility. It must be accompanied by a document outlining the problem and showing the complete contact details of the user.

The replacement of removable parts is covered for 7 years.

## VIII - GUARANTEES

This **MD** is guaranteed parts and labour against all manufacturing defects for **12 months** from the date of invoice. This guarantee does not apply to wear.

All changes or additions to the product without the express agreement of B.A. International render this guarantee null and void.

The guarantee becomes null and void if the technical instructions are not followed.

B.A. International cannot be held responsible for damage resulting from or which could result from normal wear, use, cleaning or incorrect maintenance, the non-observance of instructions for use or connection, scaling or corrosion, impurities in the water supply system or unusual chemical or electrical influences or non observance of the instructions, maintenance instructions and assembly of B.A. International and other manufacturer's instructions.

In order that guarantee requests are taken into consideration, please attach a copy of the invoice or a copy of the delivery slip to the **MD**.

EN

FR

ES

IT

DE

PO

## IX - ACCESSORIES

To be ordered from your approved distributor.

Description	Reference
Sterilisable Spraying nozzle	BA690013
Sterilisable coupling	BA690015
Needle for internal spray	BA690016
Clip for irrigation line	BA690017
Protective caps for lighting system	BA690019
B.A. International 500 ml lubricating spray	BAOIL
Spray caustic cleaner	BA690029

## X - DISPOSAL

The **MD** must be sterilised before disposal. Risk of contamination of third parties.

Comply with legislation and national standards and guidelines for disposal.

Integrated LED **MD** and adaptor: Comply with legislation, national standards and guidelines for the disposal of small electrical equipment or worn electronic equipment.

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B.A. INTERNATIONAL LTD.