Acetal Discs  
Disco de Acetal  
Disque de Acetal  
Acetal-Scheiben  
Dischi in Acetale

Acetal Schijven  
Discos de acetal  
Acetal Discs  
Acetal Discs

Distributed by (US):  
**HENRY SCHEIN INC.**  
135 DURYEA ROAD  
Melville, NY 11747 USA  
EC REP  
**HENRY SCHEIN SERVICES GmbH**  
MONZA STRASSE 2A  
D-63225 LANGEN, GERMANY  
MADE IN TRINIDAD AND TOBAGO  
Rev: 2019/09

Item Code	Product Description
9011255	ZIRLUX ACETAL MILLING DISC A1 98.5x15mm
9011260	ZIRLUX ACETAL MILLING DISC A1 98.5x20mm
9011266	ZIRLUX ACETAL MILLING DISC A1 98.5X25mm
9011256	ZIRLUX ACETAL MILLING DISC A2 98.5x15mm
9011261	ZIRLUX ACETAL MILLING DISC A2 98.5x20mm
9011267	ZIRLUX ACETAL MILLING DISC A2 98.5X25mm
9011257	ZIRLUX ACETAL MILLING DISC A3.5 98.5x15mm
9011262	ZIRLUX ACETAL MILLING DISC A3.5 98.5x20mm
9011268	ZIRLUX ACETAL MILLING DISC A3.5 98.5X25mm
9011258	ZIRLUX ACETAL MILLING DISC B1 98.5x15mm
9011263	ZIRLUX ACETAL MILLING DISC B1 98.5x20mm
9011269	ZIRLUX ACETAL MILLING DISC B1 98.5X25mm
9011259	ZIRLUX ACETAL MILLING DISC G2 (BLEACHED) 98.5x15mm
9011264	ZIRLUX ACETAL MILLING DISC G2 (BLEACHED) 98.5x20mm
9011270	ZIRLUX ACETAL MILLING DISC G2 (BLEACHED) 98.5X25mm
9012436	ZIRLUX ACETAL MILLING DISC P1 98.5 x 15mm
9012437	ZIRLUX ACETAL MILLING DISC P1 98.5 x 20mm
9012438	ZIRLUX ACETAL MILLING DISC P1 98.5 x 25mm
9012490	ZIRLUX ACETAL MILLING DISC A1 71 x 20mm
9012491	ZIRLUX ACETAL MILLING DISC B1 71 x 20mm B1
9012492	ZIRLUX ACETAL MILLING DISC G2 (BLEACHED) 71 x 20 mm

## Zirlux Acetal Instructions for Use

### SECTION 1: SYSTEM OVERVIEW

Zirlux Acetal is a pure, highly crystalline acetal copolymer resin possessing high tensile and flexural strength, fatigue resistance, and hardness. These properties, when combined with its very low moisture absorption, make the acetal material ideally suited for milling a variety of removable dental appliances, including partial denture frameworks, provisional crowns and bridges, and bite splints.

### SECTION 2: PRODUCT COMPOSITION

Ingredient	Approx % weight
Polyoxymethylene Copolymer	99 to 100
Pigments	< 1.0

### SECTION 3: TECHNICAL DATA

Property	Method	Unit	Value
Density	ASTM D792	g/cm <sup>3</sup>	1.41
Tensile Strength (New)	ASTM D638	MPa	61.00
Tensile Modulus	ASTM D638	MPa	2820.00
Flexural Strength (New)	ASTM D790	MPa	90.00
Flexural Modulus	ASTM D790	MPa	2620.00
Moisture Absorption	ASTM D570	wt%	0.22
Izod Impact, Notched	ASTM D256	ft-lb/in (J/m)	1.00

### SECTION 4: INDICATIONS

Removable prosthetic dental appliances
Partial denture clasp
Partial denture framework
Provisional crown & bridge
Bite splint

### SECTION 5: CONTRAINDICATIONS

Contraindications are all applications not included in SECTION 4 above.

### SECTION 6: PREPARATION GUIDELINES

General preparation guidelines for all removable appliances should be followed.

### SECTION 7: DESIGN OF RESTORATION

Partial denture framework (minimum thickness)

- Major connector: 0.7 mm
- Saddle area: 0.5 mm
- Clasp on anterior tooth
- 2.0 mm where clasp joins framework, thinning to 0.7 mm at tip
- Clasp on posterior tooth
- 2.5 mm where clasp joins framework, thinning to 0.7 mm at tip
- Temporary Crown & Bridge
- 0.5 mm minimum thickness
- 0.03 mm cement gap

### SECTION 8: MILLING

1. Mount the disc according to the CAM-System instructions. Set the enlargement-factor into the unit at a ratio of 1:1.
2. Start the milling unit.
3. Select appropriate milling strategy depending on the appliance to be fabricated.
4. Once the milling operation has completed, remove the disc containing the milled restoration(s).

### SECTION 9: POST MILLING

1. Remove the restoration(s) from the disc using a separating disc or suitable acrylic bur.
2. To polish the restoration, load a new polishing wheel or buff with polish suitable for thermoplastic materials and buff all areas of the appliance until all scratches are removed and high shine is achieved. Notice: All Zirlux Acetal discs are solely for use by or on the order of a dental professional. They are not for use by the general public or over-the-counter.

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