REF 958317

# LC BASE PLATES PINK

# SELEXION

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## EN Instructions for use

The light curing base plate is a light-curing resin in a pre-contoured plate for maxillae and mandibles. They are used for the fabrication of individual impressions and functional trays in crown & bridges and total prosthetic technique. They can also be used for the fabrication of bite pattern bases and dimensionally stable base plates for wax contouring. The material is ready for use and easy to handle. The material is dual cure, it will cure in day/surgery light over a long period of time however the recommended curing is via halogen or UVA light. The pleasant peppermint flavour is still apparent even after curing.

### Method of use:

**1. Fabrication of an individual tray:** Allow the light curing base plate to reach room temperature before use. Mark margins for the individual tray on the plaster model. If necessary, use a space maintainer. If space maintainer is not needed, insulate plaster model. Remove protective foil, mould plate onto the model and cut the desired rim as marked. A handle can easily be shaped from excess material. Polymerise the light curing base plates using a standard light unit (halogen or UVA light). Depending on the thickness of the tray, polymerise the reverse side for 3 - 5 min. A longer polymerisation time improves the physical properties of the material. The individual tray is immediately ready for use. A slight inhibition layer, caused by oxygen, can be removed with a tissue moistened with alcohol.


**2. Fabrication of a base plate:** Mark the margin of the planned prosthesis base on the plaster model. Block out undercuts as usual. Thinly isolate the model with an alginate isolation. Carefully place one of the light curing base plates onto the model, adjust the rims with a penknife and polymerise for 5 min in a standard halogen or UVA light unit. One single exposure is sufficient.

**Supplementary information:** For finishing the rims and the ridge, tungsten carbide metal burs are most suited. Create cross-shaped retentions on the ridge before applying wax. Then a thin layer of sticky wax can be applied.

**Storage:** Store at temperatures 4°C - 23°C. Avoid direct sunlight. Excess material can be used again if stored under light protection.

**Precautionary measures/indications:** The light curing base plate contains in an uncured condition Bis-GMA, diurethane-methacrylate, TEGDMA, BHT, peppermint oil. Do not use in cases of allergies or sensitivities. Ensure sufficient and complete polymerisation! Too short an exposure time or low intensity of light may lead to a loss of stability. Our preparations have been developed for use in dentistry. Our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice does not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for this. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

According to the EU Medical Devices Regulation, users / patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

**CB Healthcare Consulting GmbH**  
(manufacturer)  
Am Neumarkt 34  
22041 Hamburg, Germany  
Tel: +49 40 656680  
Email: [info@cbhealthcare-intl.com](mailto:info@cbhealthcare-intl.com)  
Rev.2020/07

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