# SELACRYL COLD

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Cold-curing resin for cast metal partial frames

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**EN:** SELACRYL COLD is a cold-curing resin on a methyl methacrylate basis. It is self-curing, colour stable, and exempt from tertiary amine and cadmium. SELACRYL COLD distinguishes itself for its simple processing, high mechanical properties, natural colour reproduction, and fast and secure hardening. **Preparation:** With cast metal partial frames, the wax set-up is fixed through a precast (duplicate) model in plaster, silicone or duplicating gel. The flasks, model frame and teeth are meticulously cleaned of wax residues by boiling. The plaster model must be watered to prevent air from being pressed out of the plaster into the resin during the pressure polymerization process, and the subsequent formation of air bubbles. The plaster/ acrylic resin is separated with alginate insulation. Important: the alginate insulation must be dried to avoid the white discolouration of the base area. **Mixing:** 10 g powder : 6-7 g liquid (or as desired) Pour the liquid into the mixing bowl, add the powder and stir it to a creamy and droplike consistency. Note: an excess of fluid will increase shrinkage. **Processing time:** Time to swell: 30 seconds Casting phase: 5 minutes, stable after 9 minutes. Mouldable phase: up to 14 minutes. Indicated times apply to a material and room temperature of 20-23 °C. Higher temperatures will reduce and lower temperatures will increase the indicated times. **Polymerization:** The resin mixture is poured into the flasks during the casting phase. At the beginning of the mouldable phase, the acrylic resin is solid, it no longer flows out of the flask and it can be moulded. Polymerization must be done within maximum 6-8 minutes. In a pressure vessel: approx 20-30 minutes with a pressure of 2-3 bar at 40-50 °C. For agar duplicating gel or duplicating silicone, the polymerization time is 30-40 min. **Storage and safety measures:** Material must not be used after expiration date. Storage temperature must not exceed 25 °C (77 °F). Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapors. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station. **Classification according to DIN EN ISO 20795-1 Type 2 - medical device class IIa.** According to the EU Medical Devices Regulation, users / patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred. **Cont. Methyl methacrylate.** Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. Wear protective gloves/ protective clothing/ eye protection/ face protection. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/ attention. Take off contaminated clothing and wash it before reuse. In case of fire: Use Carbon dioxide (CO<sub>2</sub>), foam or dry powder to extinguish. Store in a well-ventilated place. Keep cool. Dispose of contents/ container as hazardous waste - in accordance with local and national legislation - suitable, approved incinerator for combustible organic waste.

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
DE: Gefahr  
NL: Gevaar

Medizinprodukt  
Medical Device  
Dispositif médical  
Medisch hulpmiddel

 **CB Healthcare Consulting GmbH**  
(manufacturer)  
Am Neumarkt 34  
22041 Hamburg, Germany  
Tel: +49 40 656680  
Email: info@cbhealthcare-intl.com

 0482  
Rev.2020/01