# CROWN & BRIDGE MATERIAL 4:1

## • Instructions for use

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## INSTRUCTIONS FOR USE

### TEMPORARY CROWN AND BRIDGE MATERIAL

• Direct Application System

• Automatic mixing 2 component - bis-acrylic based provisional crown and bridge material - peroxide free and methylmethacrylate free. For production of temporary crowns, bridges, inlays, onlays, partial crowns, veneers and long-term temporary restorations.

#### Taking the impression

A study model impression should be taken using alginate or silicone before preparing for a crown or bridge or before a planned extraction. The interdental flash should be cut out to improve the stability of the temporary restoration to be fabricated at a later stage. In the case of alginate impressions, syringe the Temporary crown and bridge material as soon as possible after taking the impression.

#### Cartridge handling

The material is dispensed and mixed automatically by extruding through a mixing tip. Remove the transport seal before using a new cartridge. Then attach the mixing tip and, if required, a conturation tip. Following use, the mixing tip should be used as a seal and should be left on the cartridge until the next application. The mixing tip should only be replaced with a new one when the cartridge is to be used again.

When using for the first time, discard the material initially extruded from the mixing tip (about the size of a pea). There is an optimal mixture after initial extrusion.

#### Recommended use and elastic phase

Dry the prepared teeth and coat the preparations, surrounding tissue and any existing acrylic with a separating agent (e.g. Vaseline). Syringe Temporary crown and bridge material into the deepest section of the impression and then up to the gingival areas. Always keep the mixing tip immersed in the material to avoid bubbles. The Temporary crown and bridge material has a firm, elastic consistency approximately 2 to 3 min after mixing and can be easily removed from the patient's mouth. The curing process must be monitored intraorally (e.g. with a probe), as the temporary restoration can only be easily removed during the elastic phase.

#### Final curing and preparation

Temporary crown and bridge material cures extraorally in 3 min. The temporary restoration should be replaced in the impression during the final curing phase. Prepare and polish the temporary restoration after removing the oxygen inhibition layer with a solvent (e.g. alcohol). Commercially available eugenol-free temporary cements can be used for luting.

#### Repair of the temporary restoration

A feature of temporary restorations fabricated using Temporary crown and bridge material is its high mechanical strength. If a temporary restoration should fracture however, we recommend the following procedure:

#### Fracture of the temporary restoration shortly after fabrication:

Bond the fracture sections with freshly mixed Temporary crown and bridge ma-

terial.

#### Fracture of an older or in situ temporary restoration:

Roughen the fractured sections and prepare undercuts. Bond the prepared fracture sections with freshly mixed Temporary crown and bridge material. We recommend using a composite bond to improve the bond. Press the sections together for approx. 3 min.

#### Please note

In individual cases hypersensitivity to components of the material cannot be ruled out. Stop using the material in these cases. Avoid skin contact with the paste. On accidental contact, wash with soap and water. On contact with the eyes, rinse immediately with plenty of water and consult a doctor.

#### Technical Data

Compressive strength ≥ 200 MPa

Flexural strength ≥ 60 MPa

Diametrical tensile strength ≥ 31 MPa

Water absorption ≤ 25 µg/mm<sup>3</sup>

#### Working sequence

0 min Syringe into the impression

0 - 45 sec Insert the impression intraorally

2 - 3 min Remove from the patient's mouth

6 min Remove the oxygen inhibition layer with a solvent

6 - 7 min Preparation, contouring and polishing

#### Composition

Glass filler materials in a matrix of multifunctional methacrylates; catalyst, stabilizers, additives. Free of methyl methacrylate and peroxides. Filler content: 47 % by weight = 26% by volume. The variation width of the inorganic filler particles is between 0.02 and 2.5 µm

#### Storage

Store at room temperature (15-25 °C / 59-77 °F). Do not use after date of expiry.

#### Keep away from children! For dental use only!

According to the EU Medical Devices Regulation, users / patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

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