CYBERTECH

TEMP

CROWN & BRIDGE 4:1

Temporary Crown and Bridge Material

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Rev.2019/07

Intro-Kits

9884210 CT Crown & Bridge Material 4:1 A1, 1 cartridge 50 ml, 6 mixing tips

9884211 CT Crown & Bridge Material 4:1 A2, 1 cartridge 50 ml, 6 mixing tips

9884212 CT Crown & Bridge Material 4:1 A3, 1 cartridge 50 ml, 6 mixing tips

76 g / 50 ml

6 x

Instructions for Use

EN

TEMPORARY CROWN AND BRIDGE MATERIAL

• Direct Application System • Automatic mixing 2 component - bis-acrylic based provisional crown and bridge material - peroxide free and methylmethacrylate free. For production of temporary crowns, bridges, inlays, onlays, partial crowns, veneers and long-term temporary restorations.

Taking the impression

A study model impression should be taken using alginate or silicone before preparing for a crown or bridge or before a planned extraction. The interdental flash should be cut out to improve the stability of the temporary restoration to be fabricated at a later stage. In the case of alginate impressions, syringe the Cybertech Temporary Crown & Bridge Material as soon as possible after taking the impression.

Handling of automix syringe cartridges

Insert automix syringe into a suitable gun and remove the cap and add the appropriate 4:1 mixing tip – Blue with Orange spiral. Leave the mixing tip in situ until next use. A new mixing tip should be added when used again. On initial use there will be a small amount of wastage as the base and catalyst is extruded to product the correct mixture.

Recommended use and elastic phase

Dry the prepared teeth and coat the preparations, surrounding tissue and any existing acrylic with a separating agent (e.g. Vase-line). Syringe Cybertech Temporary Crown & Bridge Material into the deepest section of the impression and then up to the gingival areas. Always keep the mixing tip immersed in the material to avoid bubbles. The Cybertech Temporary Crown & Bridge Material has a firm, elastic consistency approximately 2 to 3 min after mixing and can be easily removed from

the patient's mouth. The curing process must be monitored intraorally (e.g. with a probe), as the temporary restoration can only be easily removed during the elastic phase.

Final curing and preparation

Cybertech Temporary Crown & Bridge Material cures extraorally in 3 min. The temporary restoration should be replaced in the impression during the final curing phase. Prepare and polish the temporary restoration after removing the oxygen inhibition layer with a solvent (e.g. alcohol). Commercially available eugenol-free temporary cements (Temporary cement) can be used for luting.

Repair of the temporary restoration

A feature of temporary restorations fabricated using Cybertech Temporary Crown & Bridge Material is its high mechanical strength. If a temporary restoration should fracture however, we recommend the following procedure:

Fracture of the temporary restoration shortly after fabrication:

Bond the fracture sections with freshly mixed Cybertech Temporary Crown & Bridge Material.

Fracture of an older or in situ temporary restoration:

Roughen the fractured sections and prepare undercuts. Bond the prepared fracture sections with freshly mixed Cybertech Temporary Crown & Bridge Material. We recommend using a composite bond to improve the bond. Press the sections together for approx. 3 min.

Please note

In individual cases hypersensitivity to components of the material cannot be ruled out. Stop using the material in these cases. Avoid skin contact with the paste. On accidental contact, wash with soap and water. On contact with the eyes, rinse immediately with plenty of water and consult a doctor.

Technical Data

Compressive strength	≥ 200 MPa
Flexural strength	≥ 60 MPa
Diametrical tensile strength	≥ 31 MPa
Water absorption	≤ 25 $\mu\text{g}/\text{mm}^3$

Working sequence

0 min	Syringe into the impression
0 – 45 sec	Insert the impression intraorally
2 – 3 min	Remove from the patient's mouth
6 min	Remove the oxygen inhibition layer with a solvent
6 – 7 min	Preparation, contouring and polishing

Composition

Glass filler materials in a matrix of multifunctional methacrylates; catalyst, stabilizers, additives. Free of methyl methacrylate and peroxides. Filler content: 47 % by weight = 26% by volume. The variation width of the inorganic filler particles is between 0.02 and 2.5 μm)

Storage

Store at room temperature (15-25 °C / 59-77 °F). Do not use after date of expiry.

Keep away from children! For dental use only!

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