C SILICONE 2

97-92688 PUTTY

97-92689 CORRECTION

97-92690 CATALYST PASTE

Instructions for use
Instrucciones de uso
Manuel d'utilisation
Gebrauchsanweisung
Istruzioni d'uso









Medical Device / Producto sanitario
Dispositif médical / Medizinprodukt / Dispositivo medicale

DEHP™ DE Healthcare Products
DEDICATED TO EXCELLENCE

 **CB Healthcare Consulting GmbH**
(manufacturer)
Am Neumarkt 34
22041 Hamburg, Germany
Tel: +49 40 656680
Rev.2019/11

C SILICONE	1-4	EN
SILICONAS POR CONDENSACIÓN	5-8	ES
SILICONES PAR CONDENSATION	9-12	FR
KONDENSATIONS VERNETZENDES SILIKON	13-16	DE
SILICONI PER CONDENSAZIONE	17-20	IT

C Silicone 2
Putty

C Silicone 2
Correction

1. ISO 4823 / ADA Spec. 19	Type 0 Putty Consistency	Type 3 Light-bodied Consistency
2. Mixing time	30"	30"
3 Total working time*	1'15	1'30
4 Time in mouth	3'15"	3'30"
5. Setting Time*	4'30"	5'00"
6. Strain in compression (min - max)	2% - 5%	2% - 8%
7. Recovery from deformation	> 98%	> 99%
8.Linear dimensional change after 24 hours	< -0,2%	< -0,7%

HYDROCOMPATIBLE CONDENSATION SILICONES (POLYSILOXANES) FOR HIGH-PRECISION IMPRESSIONS USE

C SILICONE 2 PUTTY: Polysiloxane for the first impression characterised by initial flow in the working phase and high rigidity after hardening. Recommended for the two stage technique [double impression].

C SILICONE 2 CORRECTION: Low viscosity, extremely flowable, hydrocompatible, polysiloxane. Recommended for the two stage technique [double impression].

IMPRESSION TRAY

All rigid types; DE Healthcare stainless steel impression trays are particularly recommended.

C SILICONE 2 - INSTRUCTIONS FOR USE

Dosage: Use the measuring spoon to dose the required quantity of C Silicone 2 and spread it out over your hand [note: for the correct dosage the putty should be leveled with the top of the measuring scoop]. Press the rim of the measuring scoop onto the putty once for each measure dosed. For each measure of material, apply two strips of C Silicone 2 Catalyst Paste as long as the measuring scoop [4 cm.]

Mixing: Fold and knead the putty energetically with your fingertips for about 30 sec., until the colour is even, without any stripes. Spread out the putty over the impression tray. The working time including mixing time is about 1 min. 15 sec.. Put into the mouth and wait until it is set, which will take about 3 min. 15 sec

EN

C SILICONE 2 CORRECTION - INSTRUCTIONS FOR USE

Dosage: Extrude the required quantity, observing the dosage scale on the mixing block. Then dose a quantity of C Silicone 2 equal to the length of the correction material.

Mixing: Mix energetically with the spatula, and then carefully spread the mixture out, pressing it over the mixing block to eliminate any air which may have been incorporated into the mixture. Repeat this operation until the mix is even in colour. Mixing time should ideally be about 30 sec. Apply the mixture to the impression tray with a spatula or a syringe for elastomers. Working time including mixing time is about 1 min. 30 sec. Put into the oral cavity and wait for the mix to harden, which will take about 3 min. 30 sec.

TECHNICAL DATA

Table 1 gives the characteristics of the materials:

- 1 - ISO 4823 / ADA Spec. 19
- 2 - Mixing time
- 3 - Total working time*
- 4 - Time in mouth
- 5 - Setting Time*
- 6 - Strain in compression (min - max)
- 7 - Recovery from deformation
- 8 - Linear dimensional change 24 hours

* The times for clinical use (listed in Table 1) are intended from the start of mixing at a temperature of 23°C – 73°F. Mixing at higher temperatures reduces these times. Mixing at lower temperatures increases them.

WARNING / CAUTION

- Before making the corrective impression, remove all remaining traces of oral fluids by rinsing thoroughly and drying.
- Properly close the catalyst tube after use.
- Prevent catalyst from coming into contact with skin and eyes; wash well with soap and water in the event of accidental skin contact. In the event of eye contact, rinse thoroughly with running water and seek medical attention. If substance is ingested, seek medical attention immediately.
- Please avoid contact with garments.
- In sensitive subjects, polysiloxanes may cause irritation or other allergic reactions.

EN

CLEANING AND DISINFECTING THE IMPRESSION

After rinsing well under running water, the impression can be immediately disinfected. Always follow the manufacturer's instructions.

CASTING MODELS

The impression must be cast 30 min. to 72 hours after hardening. The material is compatible with all the best gypsum available on the market. We recommend type 4 die stone and type 3 model stone.

GALVANISATION

The impressions may be galvanised with copper or silver.

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