Zirlux[®] Wax

Wax Discs

Disco in cera

Disque de cire

Wachs-Scheiben

Cera disco

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Made in Germany
Fabricado en Alemania
Fabriqué en Allemagne
Hergestellt in Deutschland
Prodotto in Germania
Gemaakt in Duitsland

HENRY SCHEIN SERVICES GMBH
MONZASTRASSE 2A
D-63225 LANGEN, GERMANY
WWW.HENRYSCH EINBRAND.COM

Distributed by (in US only)

HENRY SCHEIN INC.
135 DURYEA ROAD
Melville, NY 11747 USA

Rev. 2016/10

Size	Blue Azul Blue Blau Blu Blauw	Grey Gris Gris Grau Grigio Grijs	White, Bianco, Blanc, Weiß, Bianco, Wit
Ø98.5xH12mm	9011272	9011271	9011273
Ø98.5xH14mm	9007975-9793690	9007977-9793692	9007976-9793691
Ø98.5xH20mm	9009644-9793164	9009643-9793165	9009645-9793166
Ø98.5xH25mm	9793808	9793807	9793809
Ø100xH12mm	9011288	9011287	9011289
Ø100xH14mm	9011292	9011290	9011293
Ø100xH16mm	9011295	9011294	9011296
Ø100xH20mm	9011298	9011297	9011299

INSTRUCTIONS for dentists/dental technicians

Use: For the manufacturing of the wax scaffolding of dental implants with plastic covered crowns, through CAD/CAM systems.

1. General

The wax discs can be processed with different milling systems.

2. Storage

Protect from humidity and sunlight.

3. Area of application:

The forming wax is not designed for use in oral application. The material can be used for the preparation of dental crowns and bridges. Since it burns out leaving without residues, it is suitable for the preparation of casting moulds and ceramic moulding processes.

4. Processing

The wax milled scaffoldings can be removed with the help of appropriate cutters or cutting discs.

5a. Installation of the casting channel

The setting and the selection of the diameter of the casting channel depends on the used devices and alloys or pressed ceramics.

5b. Embedding

Line the silicon or metal muffle with a liner. Mix and embed the investment material under vacuum observing the instructions provided by the manufacturers of the investment material and alloys.

6. Disposal

Uncontaminated residues can be re-used for wax related activities. Since the product is insoluble in water, it can be disposed of together with normal household waste.

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