ZIRCONIUM DISCS / BLOCKS

• Isostatically pressed, pre-sintered

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ZIRCONIUM DISCS / BLOCKS

INSTRUCTIONS FOR USE

1. Fix the blank / block according to the CAM-System instructions. Set the enlargement-factor into the unit..
2. Start the milling unit.
3. After milling, extract the blank / block together with the contained milled restoration.
4. Remove the restoration from the blank / block with a fissure bur.
5. Remove any dust from the restoration with compressed air before sintering.
6. Position the restoration on the sinter bed of the furnace tray.
7. Program your furnace according to sintering temperature / instruction shown on the front page.
8. Run the sintering cycle.
9. Wait until cool down is completed before removing the sintered restoration.
10. Inspect the restoration in respect of flaws, wall thickness, seat and marginal fit. If necessary, make small adjustments with a water-cooled diamond bur.
11. Rinse the framework/anatomical crown with water and dry it. It is now ready for veneering. The framework is now ready for veneering /anatomical crown is ready for glazing.!

Safety precautions

- Design the restoration to maximize cross-sections.
- Carefully read the MSDS.

Indications for use

Zirconium discs / blocks are indicated for use as anatomy crown and as substructure for porcelain fused ceramic fixed dental restorations. Limitations are listed in table 1.

Table 1: Indications of use and maximum number of pontics.

Process chain	Material	Single unit crowns		Bridges, number of pontics	
		Anterior	Posterior	Anterior	Posterior
Porous blank/block is machined to enlarged framework/autonomic crown, then sintered and veneered	3Y-TZP	X	X	3	2

*X indicated; digits show maximum number of pontics. Limitation of Bridge, number of Pontics - All Zirlux-discs / blocks are solely for use by or on the order of a dental professional. They are not for use by the general public or over-the-counter.

Contraindications

Contraindications are all applications not indicated in table 1 above.

- a) Do not open the milling unit during the milling process
- b) Wear a face mask when cleaning dust with compressed air
- c) Follow your handpiece instructions carefully to avoid hand injuries
- d) Do not touch the furnace or the sintered restoration until they're completely cooled down

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