 HENRY SCHEIN®

HS-XRV SLIMLINE

- FOR DESKTOP OR WALL-MOUNTED USE

Negatoscopio Slimline

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ENGLISH

1. Preface

The HS XRV slimline is an innovative product in the x-ray film viewer/illuminator field. It can be used on desk-top or wall mounted. It features the newest LED technology, the colour temperature is more than 8,600k, this improves the light quality, making the light very even and soft. It gives better protection to the doctor's eyes, particularly if viewing the image for a long time.

2. Applications

The XRV slimline is suitable to view all sizes of x-ray film, ct film, dr film and other imaging films in hospitals, clinics, colleges and institutes. It is designed for the professional to diagnose and analyze the images.

3. Specifications

a) Size

SIZE	SCREEN SIZE L×H(mm)	OUTER SIZE L×H×W(mm)
9008520 – 36×44cm	360×440	470×545×24
9008521 – 74×44cm	740×440	850×545×24

b) Optical Specifications

- (1) No flash in the screen, luminance at centre: $>1000\text{cd/m}^2$
- (2) Optical scattering coefficient: >0.6
- (3) Even illuminating, coefficient uniformity : >0.8
- (4) Continuous working: $> 24\text{h}$
- (5) The illuminator is lit within 10 sec., the stability of luminance: < 15 minutes
- (6) Outer surface temperature after 1 hour of working : $< 60^\circ\text{c}$
- (7) Voltage: ac220v

4. Operation Instruction

- (1) The power supply is ac 220v, ensure the unit is well grounded.
- (2) Insert the plug into the outlet in the right side of illuminator.
- (3) Turn on the on/off switch at the right side, the illuminator will light immediately.

5. Warning

1. The unit should be properly earthed to prevent electrical leakage during use.
2. The plug must be well grounded, otherwise it will cause failure. If the illuminator over heats remove the plug immediately.
3. The screen is composed of an imported white polymethyl methacrylate plate, do not install the viewer in a high temperature area, as this may damage the screen.
4. Do not use the viewer in a humid place (such as a bathroom) as this is unsafe and may cause electric shocks.
5. Unplug the power plug while cleaning the viewer, or for maintenance or replacing spare parts.
6. Do not use the viewer in a flammable & explosive gas environment.
7. Remove the power plug if the unit is not used for any length of time.
8. One year warranty. Do not attempt to open for maintenance; this will negate the warranty and may cause damage.

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