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
CB Healthcare Consulting GmbH
 (manufacturer)
 Am Neumarkt 34
 22041 Hamburg, Germany
 Rev.2020/08

The CyberTech Standard system for composite restoration finishing and polishing is composed of abrasive discs and a pressure mandrel that provides agility at side inversions, as well as for disc replacements. They are offered in 9,5mm (3/8") and 12,7mm (1/2") diameters in four grits, identified by colors.

Use: Install the mandrel at the counter angle. Fit the disc onto the mandrel. In order to achieve a better result use the discs in the 4 grits, successively from the coarsest to the finest one. The product must be applied dry, with low to medium speed and with intermittent touches.

Warnings, precautions and restrictions: Professional use only. It is recommended to use 10.000 rpm for the coarsest discs and 30.000 rpm for the finest. The user will be entirely responsible for the combined use of the mandrel or the disc with similar products. To avoid an accidental swallowing of this product it's necessary the application of a total isolation technique of the operatory field. When there is proven allergy to any of its components, the product must not be used.

Connection with counter angle: The mandrel is manufactured in compliance with specifications from the rule ISO 1797.

Sterilization: The mandrel may be sterilized by autoclave at 126°C, 147kPa, 16min. Before sterilization, proceed with the cleaning using water and enzymatic detergent. The discs must be disposed after each patient. Its reuse can cause cross contamination.

Composition: Discs: Aluminum oxide, polyethylene terephthalate, synthetic rubber resin, metal eyelets, water-based pigments. **Mandrel:** Metal.

Dark green – Coarse | Light green – Medium | Yellow – Fine | Beige - X-Fine

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