External Spray E-type Airmotor Ultimate BA-603 4-hole Operation Manual

The EU directive 93/42/EEC was applied in the design and production of this medical device.

B.A. INTERNATIONAL LTD.
(manufacturer)
Unit 9, Kingsthorpe Business Centre,
Studland Road, Kingsthorpe,
Northampton NN2 6NE U.K.
Tel: +44 1604 777 700

AUTOCLAVABLE

Rev.2016/09

Indications for use

The device is intended for use in general dental applications: Drive for Straight or Contra angle handpieces for dental restoration and prophylaxis.

⚠ Caution before use


- The BA Brand Ultimate external spray airmotor is intended for use by dental professionals only.
- Read this Operation Manual carefully before using the airmotor. Be sure to follow all the operating instructions, including all cautions and warnings.
- Keep this Operation Manual for future reference.

Specifications

Back-end configuration	Applicable to 4-hole ISO standard connector
Drive Air Pressure	29-42 psi (0.2-0.29 MPa)
Rotation Speed	17,100-20,900 rpm (36 psi (0.25 MPa))
Air Consumption	less than 66 NL/min (36 psi (0.25 MPa))
Breakdown Torque	more than 1.5 N cm (36 psi (0.25 MPa))
Weight	93 grams

Note: The maximum speed may vary slightly depending on the type of tubing connected to the airmotor.

Dimensions and part descriptions:

Installation of Airmotor onto Tubing Connector

Installation

- Make sure that tubing connector is ISO standard 4-hole.
- Align the pipes at the back-end of the airmotor with the corresponding holes of the tubing, and then tighten the nut situated on the tubing by turning it clockwise.


Installation/Removal of a handpiece onto/from the Airmotor

Installation

- The BA Brand Airmotor accepts E-type handpieces with external spray. Insert the sleeve of the airmotor into the back-end of the E-type handpiece until it is securely locked accompanied by an audible snap.
- Make sure that the attachment is securely installed onto the airmotor by trying to pull it gently with the fingers before use.

Removal

- Hold the airmotor with one hand and the handpiece with the other hand, letting the hands overlap each other.
- Pull the handpiece off from the airmotor. Do not separate the overlapped hands until the lock is released.

⚠ Warning

- Never use the airmotor, if the handpiece cannot be smoothly installed onto or removed from the airmotor. In such a case, consult **BA International Ltd. or an accredited BA certified Service Centre.**

Operation

- Set the air pressure gauge to 29 psi (0.2Mpa) to 42 psi (0.29Mpa).
- Before operation, make a trial run and make sure that:

- The drive-air is not leaking from between the connecting part of the airmotor back-end and tubing connector.

• Changes in the rotating direction and rotating speed can be affected with the control ring as shown.


- When the ring is turned in "F" direction, the handpiece installed onto the airmotor rotates clockwise.

- When the ring is turned in "R" direction, the handpiece installed onto the airmotor rotates anticlockwise.

- The rotation is stopped at the centre position between "F" and "R".

- Select the rotating direction and the speed by manipulating the ring according to the treatment purposes.

External Spray System: A spray outlet is externally made. 20cm tubing with clamp for contra angle handpiece and straight handpiece installed onto the BA-603 airmotor are supplied separately as an option, if Spray, water or chip air is needed, connect the tubing to the water outlet, cutting it in an appropriate length.

Warning


- Whenever abnormality in performance, sound level, heat level, function, etc. is sensed, stop using the airmotor immediately and consult **BA International Ltd. or an accredited BA certified Service Centre.**
- Never dismantle or alter the airmotor in any way.
- Never install/remove the handpiece onto/from the airmotor whilst the airmotor is in operation.

Caution

- Never operate the airmotor with the drive-air pressure exceeding 42 psi (0.29Mpa), as the excess drive-air may cause critical damage to the internal mechanism.
- Never operate the airmotor without installing a handpiece.

Lubrication

- Use **BA ULTIMATE SPRAY.**
- The airmotor must be lubricated once a day or before autoclaving using the following procedures:
 - Remove the handpiece from airmotor.
 - Wipe off the saliva, dirt, etc. stuck to the surface of the airmotor with an alcohol-soaked soft tissue or cotton swab.
 - Keep the airmotor in a horizontal position and inject BA ULTIMATE SPRAY for three seconds through the drive-air hole.
 - Install the handpiece onto the airmotor after spray injection.
 - Make a trial run for a few seconds, moving the control ring in “F” and “R” directions several times.
 - Wipe down the excess lubricant on the surface of the airmotor with a soft tissue.

Caution

- Never use a wire brush for disinfectants for the surface cleaning.
- Never immerse the airmotor in an ultrasonic cleaner or disinfectant solution.
- Lack of lubrication after use or before autoclaving may cause premature failure of the airmotor.

Sterilization

An autoclave is recommended for sterilization.

- Before sterilization, clean the surface of the airmotor.
- Lubricate the airmotor, and then make a trial run.
- Remove the handpiece from the airmotor.
- Put the airmotor into a sterilization pouch when using a vacuum type autoclave. Do not combine it with other instruments in the same pouch.
- Autoclave the airmotor at 121 degrees C for 20 minutes or at 134 degrees C for 15 minutes.

Warning

- Never use dry heat sterilization, cold sterilization, etc.
- Never autoclave the airmotor together with other items in the same sterilization pouch.

Service and Repair

Important

This airmotor is a result of precision engineering and cannot be dismantled/repared by the third party. For service and repair, always send it to either BA International Ltd. or an accredited BA certified service centre without dismantling it.

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Autoclavable

0476

B.A. INTERNATIONAL LTD.
Unit 9, Kingsthorpe Business Centre,
Studland Road, Kingsthorpe,
Northampton NN2 6NE U.K.
Tel: +44 1604 777 700

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