BA Powersurge Premium

Instruction manual

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Manual de utilización

Manuale d'uso

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PS500 / PS500LED

BA160505 / BA160506

REF 9799547 / REF 9799546

English	3
Français	30
Español	57
Italiano	84
Deutsch	111
Português	138

Summary

1. Instructions	4	7. Adjusting parameters	16
2. Safety	4	7.1 Language selection	16
2.1 Warning symbol	4	7.2 Calibration of contra-angle in use	16
2.2 Compliant use	4	7.3 Adjustments	17
2.3 Safety instructions	5	7.4 Pedal calibration	17
3. Powersurge unit	7	7.5 Touch pad calibration	17
4. Connection / Installation	8	7.6 Restoring previous configuration / Saving current configuration	18
4.1 Connection to a local supply area	8	7.7 Software update	18
4.2 Micromotor connection / Disconnection	8	7.8 Re-initialisation of initial B.A. INTERNATIONAL configuration	18
4.3 Pedal connection / Disconnection	9	8. Using the pedal	19
4.4 Contra-angle connection	9	9. Stopping the device	19
4.5 Irrigation tube installation	10	10. Cleaning, disinfection and sterilisation	19
4.6 Pedal bar installation	11	10.1 Device maintenance	19
5. Initial use	11	10.2 Sterilisation cycle	20
5.1 Switching on	11	11. Abnormalities and errors	21
5.2 Welcome	11	12. Electromagnetic compatibility	22
5.3 Adding a user / Protocol	12	13. Guarantee	25
5.4 Changing a user	12	14. Manufacturer's responsibility	25
6. Work screen	13	15. Regulations	26
6.1 Adjustment of water flow / Speed / Torque	13	16. Description of symbols	26
6.2 Adding / Deleting / Replacing sequences	13	17. Relevant references	28
6.3 Changing sequences	14		
6.4 Graphic display	14		
6.5 Saving torque and speed data	15		
6.6 Changing the contra-angle	15		

1. Instructions

You have just purchased your new **Powersurge** dental implantology unit.

B.A. INTERNATIONAL congratulates you on your purchase and thanks you for the trust you place in it.

For complete and lasting benefit of this product, it is recommended that you read this instruction manual carefully before use.

This B.A. INTERNATIONAL product can only be used with B.A. INTERNATIONAL accessories (micromotor, pedal, calibration module, mandrel, supply cable and irrigation tubes).

The use of any other accessories is forbidden. This device is compatible with most surgical handpieces and contra-angles on the market. However, this product's technical characteristics are only valid if it is connected to an B.A. INTERNATIONAL instrument. It is an B.A. INTERNATIONAL specific connection compatible with Intramatic Lux and « E » type (standardised as per ISO 3964). All devices connected to this product must comply with CEI60601-1 : 2005 standard (if applicable).

All the precautions with  symbols must be complied with to ensure that the device is used in the best conditions possible and in complete safety.

These instructions can only be reproduced and distributed with prior agreement from B.A. INTERNATIONAL. The manufacturer reserves the right to make changes or improvements to the aforesaid equipment resulting from new technical advances. This reserve does not apply to equipment that is already installed.

2. Safety

2.1 Warning symbol

	Warning symbol.
---	-----------------

To avoid all physical and material damage, this document splits safety instructions into 3 danger levels.

	ATTENTION
	Indicates a dangerous situation which can cause material damage or mild to moderate serious injury.

	WARNING
	Indicates a dangerous situation which can cause serious / fatal injury.

	DANGER
	Indicates a great danger posed by a situation that could immediately cause serious/ fatal injury.

2.2 Compliant use

The use of this device is reserved for qualified dental health professionals, qualified or under their control. This B.A. INTERNATIONAL product is only for use in dental medicine. Any inappropriate use is forbidden.

The user is obliged :

- to only use non-defective products.
- to protect patients, third parties and themselves from all danger.
- to avoid all contamination by the product.

The user must check the safety of operation and the condition of the device before each use.

When in use, national legal requirements and the following requirements in particular, must be complied with :

- work place safety regulations.
- accident prevention regulations.

To ensure that the B.A. INTERNATIONAL product is operational at all times and to keep it in good working order, follow the maintenance, disinfection and sterilisation instructions (see section « 10 – Cleaning, disinfection and sterilisation »).

Only B.A. INTERNATIONAL technicians or technicians trained by B.A. INTERNATIONAL are authorised to repair our products..

Check settings after maintenance.

2.3 Safety instructions

2.3.1 External interference

The device is certified compliant with current electromagnetic compatibility standards.

	ATTENTION
	Users should however ensure that no electromagnetic disturbance poses an additional risk. So as not to disturb the proper operation of electromagnetic devices, it is recommended complying with the separation distances for the device (section « 12 – Electromagnetic compatibility »).

	WARNING
	The device is not designed to withstand shocks from an electrical defibrillator.

2.3.2 Electrical connection

	DANGER (ELECTRIC SHOCK / ELECTROCUTION)
	Non-compliant electrical connection Check that the supply voltage is the same as that shown on the control unit. Check the supply cable before each use. Never stick or remove a plug from a wall socket with damp hands. Do not connect USB peripheral having an external power supply. Any other use or connection of devices to USB port may create new risks and may require additional risk analysis.

2.3.3 Environment

	ATTENTION
	Do not tilt the device. Ensure that cables do not prevent free movement of individuals. Do not use the device near ionising radiation. The device must be connected to a nearby wall socket which must be easily accessible. The required local supply cables include a socket supporting at least 125% of the nominal voltage for the equipment. Use it on a maximum impedance network of 4 Ω.

	DANGER Do not use this product in facilities that are at risk of explosion. Do not expose the device to water projections or water spray. Do not insert metal items in the device in order to avoid the risk of electrical shock, short-circuit or the emission of dangerous substances.
	DANGER (ELECTROCUTION) Accidental leak of fluid into the box. Check and ensure the water tightness of containers and pipes before each use. Should fluid leak into the device, do not touch it and immediately disconnect the electricity supply cable. Ensure that the surface of the device is completely dry before re-connecting the electricity cable.

2.3.4 Using the device

	ATTENTION Do not move the device when in use. We strongly recommend that torque is restricted so as not to exceed the manufacturer's recommendations for the instruments. This is to reduce the risk of injury and instruments being damaged. When the micromotor rotation direction is changed, this is recorded.
	WARNING Switch the device off as soon as the treatment is finished. Disconnect it from the electricity supply when not in use for a long time.
	ATTENTION There are mobile parts in the peristaltic pump. Do not insert fingers inside the pump cabinet.

2.3.5 Storage conditions

	ATTENTION Recommended conditions : keep away from humidity and in ambient temperature. Keep away from the sun, sparks and all sources of ignition. Temperature : from -20°C to 70°C. Air humidity levels : from 5 % to 95 %. Incompatible materials : acids and other solvents.
---	--

2.3.6 Malfunction

In the event of malfunction (screen that glitters, noisy pump motor, vibration in the micromotor) :

- Immediately stop using the device.
- Try to locate or eliminate the problem using the description (section « 11 – Abnormalities and errors »).

→ If it is not possible to locate or eliminate the problem using the description, switch the device off and call your B.A. INTERNATIONAL approved repair center.


	DANGER
	Never try to repair it alone.

3. Powersurge unit

Fig. A

1. Touch screen.
2. Micromotor connection socket.
3. Peristaltic pump cabinet.
4. Integrated calibration module.
5. Micromotor.
6. Micromotor connection cable.
7. Motor support.
8. Pedal.
9. Supply switch (On/ Off).
10. Supply socket.
11. USB socket for storage peripherals (USB key).
12. USB socket for After-Sales Dept. use.
13. Pedal socket.
14. Calibration mandrel.

There are 2 versions of Powersurge units available :

- One version without light.
- One version with an integrated LED in the contra-angle (see section « 17- Relevant references »).

This product's software is a « multi-user » and/or « multi-protocol of implant », i.e. it is designed so that several users can use the same resources. We can create up to 10 users or protocols that can hold a maximum of 10 sequences.

Technical characteristics

- Input supply : from 100 to 240 V 50-60Hz.
- Absorbed power : 200 W.
- Operational mode : intermittent S3 10 % (1 min ON/ 9 min OFF).
- Fusibles Ø 5x20 T3, 15AH 250 V.
- Motor regime : 300 to 40 000 rpm.
- Motor torque : 7 N.cm (80 N.cm max with a **Powersurge** PS20X, L, LED - 20 :1 contra-angle).
- Pump : 5-wheel peristaltic.
- Flow : 120 mL/min maxi (with B.A. INTERNATIONAL irrigation lines).
- Motor : BF applied part.
- Coupling type : B.A. INTERNATIONAL specific Intramatic Lux compatible and « E » Type.
- IP2x box .
- IPx8 pedal.
- Conditions of use : temperature: +10°C to +40°C/ Relative humidity: 15 to 80 %.

Units displayed and meanings.

N.cm = torque [Newton per centimeter].

Rpm = speed [Turns per minute].

4. Connection / Installation

One must be familiar with all the safety instructions before installing and using the device [see section « 2 – Safety »].

4.1 Connection to a local


Only use the supply cable supplied by the manufacturer.


Grounding reliability can only be achieved when power supply is connected to an equivalent receptacle marked « Hospital Only » or « Hospital Grade ».

Check that the supply voltage is the same as that shown on the control unit and that the supply switch is in the 0 position.

4.2 Micromotor connection/ Disconnection

Line up the red dot on the micromotor with the red dot on the connection on the front side of the panel and clip :



To disconnect it, remove the safety catch as follows :
Pull the mobile part of the connection in the direction of the arrow.


4.3 Pedal connection / Disconnection

To connect the pedal :

Have the pedal socket as shown
(locking strip towards the bottom) and clip.

To disconnect the pedal :

Push the strip in the direction of the arrow
to disconnect the socket.


4.4 Connection of contra-angle

Clip the contra-angle to the motor as follows :

For proper connection, line up 2 by 2 the following items as follows.
[a] LED connector, **[b]** light claw.

Reverse the procedure to disconnect the contra-angle :

4.5 Installation of the irrigation line

Towards the
contra-angle


Towards
the bag


1. Press the « PUSH » logo to open the box cabinet.
2. Install the irrigation line by following the IN/OUT directional arrows shown and place the 2 plastic parts (with larger diameter) outside of the cabinet.
3. Close the cabinet.
4. Connect the screw nozzle (LUER) to the fluid bag or perorate the cap of the bottle. Suspend it from the bracket for this purpose.

Guide the irrigation line along the micromotor cable using the attachment clips as follows

4.6 Pedal bar installation

1. Remove the two screws having discs.
2. Insert the bar in accordance with the illustration.
3. Replace the screws + discs and screw moderately.

5. Initial use

One must be familiar with all the safety instructions before installing and using the device [see section E « 2 – Safety »].


5.1 Switching on

Put the supply switch in the position (**Fig.A, rep.9**).

5.2 Welcome

The welcome screen allows the selection of the user or the wished protocol.
This is the description of the welcome screen according to the representation **Fig.1** :

Fig. 1

1. Users or Protocols.
2. Button and name of selected user/ protocol.
3. User button / Inactive protocol.
4. Access settings.
5. Validation button for user selection.
6. Software version.

If the LCD screen does not come on, see section « 11 – Abnormalities and errors ».


5.3 Adding a user / protocol

 It is possible to configure up to 10 users / protocols that can include up to 10 sequences.

Example : we are going to add the « DOCTOR » user associated with the icon 

	Press the « settings » button (Fig.1, rep.4 or Fig.3, rep.11).
	Press the button (Fig.9, rep.1), to access the users/ protocols version screen.
	Press the « add a user » button (Fig.2 rep.5).
	Select an icon for the new user/ protocol and validate 
	Name the new user/ protocol and validate (12 letters maxi.).

Fig. 2

Here is the description of the screen « modification of users/ protocols » according to representation Fig.2.

- 1. Edition of sequences of user / active protocol.
- 2. Modification of the icon (double-click).
- 3. Delete user / protocol.
- 4. Adding a user / protocol.
- 5. Modification of user's name / active protocol.

5.4 Changing a user

	Press the icon for the active user (Fig.3, rep.1) to access the welcome screen (Fig.1) where the list of users already programmed is shown.
	Select the user of your choice and validate 

 If the user of your choice is not in the list, select the settings icon  (Fig.1, rep.4 or Fig.3, rep.11) and refer to sections « 5.3 – Adding a user/ protocol »

6. Work screen

Select a user protocol and validate (**Fig.1, rep.5**) to access the work screen (**Fig.3**).


When switching on the device, the welcome screen only stays active for 10 seconds before toggling to the work screen (**Fig.3**) associated with the active user / protocol.

When functioning, only the main data are displayed on the screen : Flow / Torque / Speed.

This is the description of the « work screen » according to the representation **Fig.3**.

Fig. 3

1. Active User / Protocol and access to welcome screen.
2. Reduction ratio for active contra-angle and choice of it.
3. Activate or deactivate the brightness of the contra-angle.
4. Instantaneous speed of the tool if the motor is active or speed limit when stopped.
5. Instantaneous torque if the motor is active or maxi. reached when the motor is stopped.
6. Active sequence of all the sequences included in the active protocol.
7. Programmed torque limit.
8. Access to graphical display.
9. Edition of sequences of the current protocol.
10. Rotation direction of the motor and reversed direction.
11. Access settings.
12. Adjust the flow of the peristaltic pump.
13. Irrigation level.

Tactile zone 13, 4, 5, 12, 7 : Access settings Flow/ Speed/ Torque.

Press the centre of the screen (**Fig.3**) to have access to the settings of the 3 main data: irrigation flow, maximum speed and maximum torque.


With the work screen, it is also possible to adjust the flow of peristaltic pump with a short press on button A of the pedal (see section « 8 – Using the pedal »).


6.1 Settings Flow / Speed / Torque

This is the description of the screen « settings Flow/ Speed / Torque » according to the representation **Fig.4** :

Fig. 4

1. Flow of peristaltic pump.
2. Decrease the value.
3. Adjusted flow.
4. Increase the value.
5. Setting cursor.
6. Speed of micromotor.
7. Set maxi. speed.
8. Torque limitation of the micromotor.
9. Maximum reachable torque.
10. Cancel.
11. Validate.

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6.2 Adding / Deleting / Replacing sequences

Press the icon (**Fig.2, rep.1 or Fig.3, rep.9**) to access the sequence management screen (**Fig.5**).

6.5 Saving torque and speed data

The procedure to follow for saving « torque and speed data » is as follows :

	Press the button (Fig.3, rep.8) .
	Press the button (Fig.6, rep.7) .
	Name the file and then validate.

Data « torques and speeds » are located in the index "ImplMotCurve".

A file called « name » is created during backup. « name » is corresponding to what you captured.

A file called « name.csv » is created there during backup, if the file already exists, an incremental file is added : "Name001.csv".

It is a text file format that can be imported into a spreadsheet (Microsoft Excel for example).

Those files can be directly imported into a spreadsheet.

Fig. 7

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The following data are measured every 0.25 s :


TIME : display time (in ms) that passed since the device is switched on.

TORQUE : display torque in 1/100 of N.cm (or in 1/100 of % if the reduction ratio is different from 20:1).

SPEED : display motor speed, in rpm.

.png images format **(Fig.8)** are also created in the index « name », one file for each sequence where recording of torque and speed **(Fig.5, Rep.5)** is activated.

Fig. 8

6.6 Changing a contra-angle

This is the list of available contra-angles :

1:3

1:1

5:1

16:1

20:1

30:1

256:1

Press the icon for the reduction ratio for the active contra-angle **(Fig.3, rep.2)**.

Note that the speed ranges at the exit of the tool are displayed to the side of the selected contra-angle.

7. Adjusting parameters

This is the description of the screen « configuration » according to the description **Fig.9** displayed after pressing the button .

Fig. 9

- 1. Users or Protocols.
- 2. Language selection.
- 3. Calibration of used contra-angle.
- 4. Adjustment of LED lighting.
- 5. Adjustment of sound level.
- 6. Advanced settings.
- 7. Pedal calibration.
- 8. Touch pad calibration.
- 9. Save/ transfer of the parameters.
- 10. Software update.
- 11. Re-initialisation of initial configuration.
- 12. Return to previous screen.

7.1 Language selection

	Press the « settings » button (Fig.1, rep.4 or Fig.3, rep.11).
	The button (Fig.9, rep.2) accesses the language selection screen.
	Then select the required flag and validate.

7.2 Calibration of contra-angle in use

 To make the control and display of the torque value more accurate, the contra-angle needs to be calibrated.

This provides the motor with information about output and efficiency of the contra-angle in use.

The test is carried out in 2 main phases :

- First, the contra-angle is made to rotate with no load.
- The contra-angle must then be connected to the specific calibration module. This module simulates an average intervention in the mouth. The measurement is therefore made in conditions as close as possible to real operation.
- Before using the calibration module, it MUST be lubricated using the B.A. INTERNATIONAL spray Ref. **BAOIL** equipped with efitting nozzle.
- Place the adaptor over the connection orifice.
- Maintain the spray adaptor in contact on the calibration module.
- Spray for 1 second.

 The calibration module is not sterilizable.

 As the capacities of the contra-angles can vary significantly, calibration of the contra-angle is recommended after maintenance and sterilization, with assistance of a second person to bring the gauge to a standstill (non sterile) and retrieve the mandrel after calibration.

	Press the « settings » button (Fig.1, rep.4 or Fig.3, rep.11).
	Press the button Fig.9, rep.3 .

	Connect the 1/20 contra-angle [equipped with the specific connection tool] to the motor. Important : do not connect the calibration module at this stage. The motor must first perform a no-load test. When the contra-angle is connected, press the pedal lever. The no-load test is run. The motor increases in speed for approximately 10 seconds. At the end of this first test, a screen is displayed asking you to connect the contra-angle load module. When connected, validate again using the pedal lever. The load test is run correctly. Maintain the tool calibration module connection correctly. The motor increases in speed for approximately 10 seconds. When it is terminated, the motor stops automatically and the screen displays « contra-angle calibrated ». You are then automatically returned to the main parameterization screen. Your contra-angle is now calibrated.
---	--

 After being used 100 times or after 2 years of running, the calibration device has to be returned to the B.A. INTERNATIONAL After-Sales Service. It will then be checked and restandardized by B.A. INTERNATIONAL.

7.3 Adjustments

It is possible to adjust the light intensity of the contra-angle, the sound, to manage the alarms and the temporizations with the following procedure :

	Press the button « <i>adjusting parameters</i> » (Fig.1, rep.4 or Fig.3, rep.11).
	Press on the button of your choice (Fig.9, rep.4 or Fig.9,rep.5 or Fig.9,rep.6).
	Press on the buttons  and/ or  move the cursor on the bargraph.
	Press on the button .

7.4 Pedal calibration

	Press on the button « <i>settings parameters</i> » (Fig.1, rep.4) or (Fig.3, rep.11).
	Press on the button (Fig.9, rep.7).
	Keep on pressing on the pedal's lever and press on the button « <i>validate</i> » 
When it is terminated, a screen displays « <i>calibration successful</i> »	

7.5 Touch pad calibration

If the touch pad is not working, it may need to be adjusted.

To do this :

	Press the « <i>settings</i> » button (Fig.1, rep.4) or (Fig.3, rep.11).
	Press the button (Fig.9, rep.8) to access the touch pad calibration screen.
	Press on the indicated areas.
When it is terminated, a screen displays « <i>calibration successful</i> »	

7.6 Restoring previous configuration / Saving current configuration

This function saves all the customisation done in order to be able to subsequently restore it.

The following can then be saved :

- Various users, protocols and sequences created with the associated parameters.
- Various adjustments made.

Connect the USB key to the rear side of the box.

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Connect the USB key to the rear side of the box.

	Press the « <i>settings</i> » button (Fig.1, rep.4) or (Fig.3, rep.11).
	Press the button (Fig.9, rep.10).
	Validate to return to the « <i>settings</i> » screen.


Updating the motor does not result in the deletion of any previously configured parameters.

7.8 B.A. INTERNATIONAL initial configuration

	Press the « <i>settings</i> » button (Fig.1, rep.4) ou (Fig.3, rep.11).
	Press the button (Fig.9, rep.11).
	Validate to return to the « <i>settings</i> » screen.

8. Using the pedal

The lever gradually adjusts the speed according to the motor during the operational phase.

A short press on the button « A » is a 20% increase in the speed of the flow of the peristaltic pump.

A long press on the button « A » starts the pump at the programmed flow, without using the micromotor, or the tool.

A short press on the button « B » inverts the rotation direction of the tool.

A long press on the button « B » increases the limitation of the torque in increments of 5 N.cm.

A short press on the button « C » goes on to the next sequence for the active protocol.

A long press on the button « C » goes to the previous sequence for the active protocol.

If the last protocol sequence is active, **a short press on button « C »** makes the first sequence of this protocol active.

The first protocol sequence is active, **a long press on button « C »** makes the last sequence of this protocol active.

9. Stopping the device

At the end of the dental surgery :

- Put the supply switch in the 0 position.
- Remove the tool from the contra-angle.
- Remove the bottle or bag from the bracket.
- Remove the perforator of the irrigation line from the bottle or bag of physiological serum or sterile water.
- Remove the clips for irrigation lines.
- Disconnect the contra-angle from the micromotor.

10. Cleaning, disinfection and sterilisation

10.1 Device maintenance

Single use	Non sterilisable	Sterilisable in an autoclave at 134°C for 18 minutes	
		Micromotor sterilisation cap	
		Sterilisable non immersible (= disinfection)	Sterilisable and immersible (+ ultrasound)
→ The bags of physiological fluid → The irrigation lines → Screen protection → Micromotor protective sheath	→ Control unit → Pedal → Calibration module	→ Micromotor + cable	→ Bracket → Motor support → Mandrel → Fixing clips for the irrigation line

*All of this information is summarised in a maintenance poster supplied with the motor.

The sterilisable products are supplied non sterile and need to be decontaminated and sterilised prior to first use and immediately after each use. Use disinfectant products that comply with the instructions given by the manufacturer.

B.A. INTERNATIONAL sterilisable medical products can withstand a maximum temperature of 138°C.

Remove the irrigation line and place it in a container for soiled medical devices (single use irrigation line).

Do not use the spray directly on the device to clean it; use disinfectant wipes instead.

10.2 Sterilisation cycle

11. Abnormalities and errors

This section is intended to help users should a problem occur.

In the event of abnormal operation (screen that glisters, noisy pump motor, vibration in the micromotor), immediately stop using the device and call your B.A. INTERNATIONAL approved repair center. If it is not possible to locate or eliminate the problem using this description, or if the disfunction symbols remain active, switch off the device and call your B.A. INTERNATIONAL approved repair center. Under no circumstances try to dismantle the control unit, motor or control pedal for the **Powersurge** unit.

Never try to repair it alone.

Abnormality screen	Analysis of disfunction
	The pedal is not detected. → Check that the pedal is properly connected to the box and check the status of the connector. → Press on the screen to make the safety signal go away and to go back to the previous or initial screen.
	The motor is not connected. → Check that the motor is properly connected to the box and check the status of the connector. → Press on the screen to make the safety signal go away and to go back to the previous or initial screen.
	During the calibration process, the loading module is not connected to the contra-angle. → Connect the loading module and validate to restart the operation. → Press on the screen to make the safety signal go away and go back to the previous screen.
	Excess speed of the motor detected. → Press on the screen to make the safety signal go away and go back to the working screen.

Dysfunction	Causes	Solutions
The screen does not come on.	→ The device is switched off. → The fuses are burnt. → The wall supply socket is defective. → The supply cable is defective.	→ Check that the device is properly connected to the supply and that the On/ Off button is in position 1. → Replace the fuses. → Move the device. → Replace the cable.

Dysfunction	Causes	Solutions
The tool is not turning.	→ Incorrect insertion of the tool. → Broken contra-angle. → Blocked motor.	→ Check the installation of the tool. → Change the contra-angle. → Check the status of the motor.
No spray or low flow or water tightness problem.	→ The pump flow is 0% or is too weak. → Incorrect installation of irrigation line. → Irrigation line defective or jammed.	→ Adjust the flow of the pump from the working screen or the pedal. → Check the positioning of the line. Check that the run of direction matches the IN/ OUT indications shown on the pump cabinet. → Replace the irrigation line.
A red screen appears indicating an error.	→ The data returned by the motor is incorrect. → There may be several reasons for this : - contra-angle malfunction - micromotor malfunction - motor malfunction	→ Switch on the motor again. → If the problem persists, note the error number and return the contra-angle unit, micromotor and motor to your after-sales service centre.

12. Electromagnetic compatibility

This product requires specific precautions with respect to electromagnetic compatibility. It must be installed and fitted according to the instructions in the section « 4 – Connection / Installation ».

Some types of mobile telecommunications devices such as mobile phones can interfere with this product.

The separation distances recommended in this section must be followed.

This product must not be used near or on any other device. If this is unavoidable, one must check its proper operation in the usage conditions prior to use.

The use of accessories other than those specified or sold by B.A. INTERNATIONAL as replacement parts, can result in increased emissions or a reduction in the immunity for this product.

Directives and manufacturer's declaration – Electromagnetic emissions

The **Powersurge** unit is intended to be used in the electromagnetic environment specified below. The customer or user of the **Powersurge** unit must ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - directives
RF CISPR 11 emissions	Group 1	The Powersurge unit uses RF energy only for its internal functions. Consequently, its RF emissions are very low and are not likely to cause interference to an electronic machine nearby.

Emissions test	Compliance	Electromagnetic environment - directives
RF CISPR 11 emissions	Class B	The Powersurge unit is suitable for use in all facilities, including domestic facilities and those directly connected to the low voltage public electricity network supplying buildings for domestic use.
Harmonic emissions CEI61000-3-2	Class A	
Voltage fluctuation / Flicker twinkling CEI61000-3-3	Compliance	

Directives and manufacturer's declaration – Electromagnetic immunity

The **Powersurge** unit is intended to be used in the electromagnetic environment specified below.
The customer or user of the **Powersurge** unit must ensure that it is used in such an environment.

Immunity test	Test level CEI 60601	Compliance level	Electromagnetic environment - directives
Electrostatic discharges CEI 61000-4-2	± 6 kV on contact ± 8 kV in the air	± 6kV on contact ± 8kV in the air	Floors must be wooden, concrete or ceramic tiles. If floors are covered with synthetic materials, the relative humidity must be at least 30 %.
Rapid transition into salvos CEI 61000-4-4	± 2 kV for electricity supply cables ± 1 kV for input/output cables	± 2 kV for electricity supply cables ± 1 kV for input/output cables	The quality of the electricity supply network must be that of a typical commercial or hospital environment.
Transitory excess voltage CEI 61000-4-5	± 1 kV between phases ± 2 kV between phases and earth	± 1 kV between phases ± 2 kV between phases and earth	The quality of the electricity supply network must be that of a typical commercial or hospital environment.
Brief cut in voltage and variations in voltage in input electricity supply cables CEI 61000-4-11	<5 % UT (>95 % cavity of UT) for 0,5 cycle 40 % UT (60 % cavity of UT) for 5 cycles 70 % UT (30 % cavity of UT) for 25 cycles <5 % UT (>95 % cavity of UT) for 5 s	<5 % UT (>95 % cavity of UT) for 0,5 cycle 40 % UT (60 % cavity of UT) for 5 cycles 70 % UT (30 % cavity of UT) for 25 cycles <5 % UT (>95 % cavity of UT) for 5 s	The quality of the electricity supply network must be that of a typical commercial or hospital environment. If the user of the Powersurge unit needs continuous operation during cuts in supply, it is recommended supplying the Powersurge unit from energy supply with no cut-out or with a battery.
Magnetic field for the frequency of the electricity network [50/60Hz] CEI 61000-4-8	3 A/m	3A/m	The magnetic fields at the frequency of the electricity network must have the levels characteristic of a representative place located in a typical commercial or hospital environment.

NOTE : U_i is the alternative voltage before applying the test level

Directives and manufacturer's declaration – Electromagnetic immunity

The **Powersurge** unit is intended to be used in the electromagnetic environment specified below.
The customer or user of the **Powersurge** unit must ensure that it is used in such an environment

Immunity test	Test level CEI 60601	Compliance level	Electromagnetic environment - directives
RF disturbances CEI 61000-4-6	3 Veff from 150 kHz to 80 MHz	3 V	Portable and mobile RF communications devices must not be used near any part of the Powersurge unit, including the cables, for the recommended separation distance, calculated from the equation applicable to the frequency of the emitter. Recommended separation distance $d=1,2 \sqrt{P}$
Radiated RF disturbances CEI 61000-4-3	3 V/m from 80 MHz to 2,5 GHz	10 V/m	$d=0,4\sqrt{P}$ from 80 MHz to 800 MHz $d=0,7\sqrt{P}$ from 800 MHz to 2,5 GHz where P is the power characteristic for maximum output of the emitter in watts (W), according to the manufacturer of the emitter and d is the recommended separation distance in metres (m). The strengths of the fixed RF, determined by an electromagnetic investigation on site a, must be less than the compliance level, in each frequency range b. Interference can occur near the device marked with the following symbol: 

NOTE 1 : to 80 MHz and at 800 MHz, the highest frequency range applies.

NOTE 2 : these directives do not apply in all situations. Electromagnetic propagation is affected by absorption and by the reflections of structures, objects and individuals.

^a The strengths of fixed emitter fields, such as base stations for radiotelephones (cellular/cordless) and terrestrial mobile radios, amateur radio, AM and FM broadcasting, and TV broadcasting, cannot be theoretically tested with accuracy. To assess the electromagnetic environment due to fixed RF emitters, an electromagnetic investigation must be considered on site. If the strength of the field, measured in the location where the **Powersurge** unit is used, exceeds the RF compliance level applicable above, the **Powersurge** unit should be observed to check that it is working normally. If abnormal performance is found, additional steps may be necessary such as redirecting or repositioning the **Powersurge** unit.

^b For the frequency range from 150 kHz to 80 MHz, the strengths of the field must be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications devices and the Powersurge

The **Powersurge** unit is intended to be used in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **Powersurge** unit can help to prevent electromagnetic interference by keeping a minimum distance between the portable and mobile RF communications device (emitters) and the **Powersurge** unit, as recommended below, according to the maximum emission power of the communications device.

Maximum output power assigned for the emitter W	Separation distance according to the frequency of the emitter (m)		
	From 150 kHz to 80 MHz $d=1,2/\sqrt{P}$	From 80 MHz to 800 MHz $d=0,4/\sqrt{P}$	From 800 MHz to 2,5 GHz $d=0,7/\sqrt{P}$
0,01	0,12	0,04	0,07
0,1	0,38	0,13	0,22
1	1,2	0,4	0,7
10	3,8	1,3	2,2
100	12	4	7

For emitters whose maximum emission power assigned is not given above, the recommended separation distance d in meters (m) can be estimated by using the equation applicable to the frequency of the emitter, where P is the characteristic of the power of maximum emission of the emitter in watts (W), according to the manufacturer of the latter.

NOTE 1 : to 80 MHz and at 800 MHz, the separation distance for the highest frequency range applies.

NOTE 2 : these directives do not apply in all situations. Electromagnetic propagation is affected by absorption and by the reflections of structures, objects and individuals.

13. Garantie

12 month GUARANTEE.

- This device is guaranteed parts and labour against all manufacturing defects for 12 months from the date of invoice.
- This guarantee does not apply to wear and tear parts.
- All changes or additions to the product without the express agreement of B.A. INTERNATIONAL render this guarantee null and void.
- The guarantee becomes null and void if the technical instructions are not followed.
- B.A. INTERNATIONAL cannot be held responsible for damage resulting from or which could result from normal wear, use, cleaning or incorrect maintenance, the non-observance of instructions for use or connection, scaling or corrosion, impurities in the water supply system or unusual chemical or electrical influences or non observance of the instructions, maintenance instructions and assembly of B.A. INTERNATIONAL and other manufacturer's instructions.
- In order that guarantee requests are taken into consideration, please attach a copy of the invoice or a copy of the delivery slip to the device.

14. Manufacturer's responsibility

This leaflet does not substitute leaflets for peripheral products.
The manufacturer reserves the right to change the device and/or the user manual without warning.

The following points, if they are not followed, relieve B.A. INTERNATIONALE of the responsibility :

- All the manufacturer's recommendations must be followed during installation (network voltage, electromagnetic environment...).
- The device must be used in an electrical facility compliant with current regulations.

- The use of this device must stay in the context of the operation for which it was designed.
- No changes, extensions and settings must be made by third parties or any other individual.
- Opening and repairs can only be carried out by individuals trained by B.A. INTERNATIONAL and who have been accredited.
- Only B.A. INTERNATIONAL original replacement parts must be used.
- The product must be used in accordance with the assembly, service and maintenance instructions..
- All the instructions contained in this document must be followed.

15. Regulations

According to ISO11498 and ISO7494 standards, the device is in the category of dental treatment devices.

This medical device is classed IIa as per the European directive 93/42/CE.

This equipment is manufactured in compliance with the following current standard IEC60601-1.

This equipment was designed and manufactured according to ISO 13485 quality assurance system.

As an electric and electronic equipment, the disposal of the product must be done by a specialised collection facility for the removal and recycling or destruction (in particular in the European market, with reference to directive 2002/96/CE dated January 23rd, 2003).

When the product becomes obsolete, we recommend you contact your supplier of dental equipment to determine the action to take.


16. Description of symbols

	Alternative voltage.
	BF type device.
	Attention : follow the instructions for use.
	Class II device.
	Manufacturer.
	Class I or IIa medical device compliant with European directive 93/42CEE.
	Gost certification.
	Thermoisinfectable device.
	Sterilisation in an autoclave except for packaging.
	Do not use if package is damaged.
	Protect from light.
	Authorised temperature limit.
	Fragile.

	Keep dry.
	Handle with care.
	Device expiry date.
LOT	Device manufacturing batch number.
SN	Serial number.
50 kg MAX 	Authorised stacking load 50 Kg.
	Humidity limit.
	Atmospheric pressure limit.
REF	Commercial reference number for the device.
	Transport upright ; with the top in the direction of the arrow.
STERILE	Sterilised device.
	Do not re-use, single-use device.
	Do not throw in domestic rubbish.
	Do not move.
IP2x	IP2x Protected from solid foreign bodies which Ø is 12,5 mm and more.
IPx8	IPx8 Protected from the prolonged effects of water.
	No sterilization.
	Desinfecting wipe cleaning.
	Soaking.
	Ultrasonic.
	Rinsing.

	Drying.
	Sterilization during 18 minutes.

17. Relevant references

Accessories	Reference	Description
Pack of 10 irrigation lines for bag or bottle*.	BA690010	
10 attachment clips for tubing.	BA690008	
Micromotor sterilisation cap.	BA690004	
Calibration module box.	BA690000	
Mandrel for calibration module.	BA690003	
Motor support.	BA690009	
Bracket.	BA690006	

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BA Powersurge Premium

Accessories		Reference	Description
Powersurge Implantology 20 :1 LED contra-angle <i>Without depth stop system.</i>	☀	BA660571	
Powersurge Implantology 20 :1 contra-angle <i>Without depth stop system.</i>		BA660070	
Powersurge Surgery LED contra-angle.	☀	BA660572	
Powersurge Surgery contra-angle.		BA660071	
Powersurge Surgery LED handpiece.	☀	BA660574	
Powersurge Surgery handpiece.		BA660072	
Powersurge contra-angle – latch type head.		BA660073	
Powersurge contra-angle with push button head.		BA660074	
Micromotor without light (with BA690001 cable).		BA640070	
LED micromotor (with BA690001 cable).	☀	BA640571	
Micromotor connection cable.		BA690001 (2 m cable)	-

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B.A. INTERNATIONAL LTD.