



## Instruction for Use

### **Absorbent Paper Points**

Sterilized endodontists points used in root canal therapy to remove moisture from the root canal or to introduce medicaments.

### **Composition**

100% Paper

### **Application**

Remove moisture from the root canal or Introduce medicaments.

### **Precaution**

- Professional use only.
- Do not autoclave.
- Single-use device / Re-use is not allowed.

### **Shelf Life**

- 4 years from the date of manufacture

### **Symbols EN ISO 15223-1:2016**



Use by Date



Lot number



Manufacturer



Single Use



Do not use if  
package is  
damaged



Sterilized through  
irradiation



Medical Device

### **Warranty**

CB Healthcare Consulting GmbH will replace the product that is proven to be defective.

CB Healthcare Consulting GmbH does not accept liability for any damage or loss, direct or consequential, stemming from the use of or inability to use the product described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use.

The user assumes all risk and liability in connection therewith.